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Stand: 10/2022  
(Rev. 2.0)

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Status: 10/2022  
(Rev. 2.0)

## Table of contents

1	Foreword.....	4
2	General information.....	4
2.1	Symbols used .....	4
2.2	Type plate.....	6
2.3	Proof of standards.....	7
3	Safety instructions.....	8
4	General product description .....	9
4.1	Intended use .....	9
4.2	Indication .....	9
4.3	Contraindication.....	9
4.4	Equipment features.....	10
5	Assembly information .....	11
5.1	Basic information on installation .....	11
5.2	Dimensional sketch of the movement space .....	11
5.3	Dismantling the care bed .....	12
6	Operation .....	12
6.1	Operating the grab rails .....	12
6.2	Shut-off function .....	12
6.3	Operating the functions .....	13
6.4	Operating the Trendelenburg function .....	13
6.5	Operating the rotary function .....	14
6.6	Operating the brake rollers .....	15
6.7	Patient lifting pole with grab handle.....	15
6.8	Operating instructions.....	16
6.9	Mattresses approved for use .....	17
7	Ambient conditions .....	17
7.1	Storage conditions.....	18
7.2	Operating conditions .....	18
8	Technical data subject to change!.....	18
9	Materials used.....	19
10	Service and care .....	19
11	Useful life of the product .....	19
12	Disinfection .....	19
12.1	Specifications of the cleaning agents and disinfectants:.....	20

13	Malfunctions and their rectification.....	20
14	Recommended accessories.....	21
15	Maintenance .....	21
15.1	Legal basis .....	21
15.2	Maintenance intervals.....	21
15.3	Spare parts .....	21
15.4	Notes on documentation .....	22
16	Re-use .....	22
17	Waste disposal .....	22
17.1	Disposal of the appliance .....	22
17.2	Disposal of electrical components.....	22
17.3	Disposal of packaging.....	22
18	EC Declaration of Conformity .....	24

## List of illustrations

Figure 1:	Exemplary type plate .....	6
Figure 2:	Dimensional sketch of the movement area of the Rotadorm Mobile .....	11
Figure 3:	Removable patient handle.....	12
Figure 4:	Shut-off box .....	12
Figure 5:	Patient hand control.....	13
Figure 6:	Bed in the standing position.....	14
Figure 7:	Castors in braked and unbraked condition .....	15
Figure 8:	Gallows mount at the head end.....	15
Figure 9:	Adjustable handle .....	16

## List of tables

Table 1:	Verification of standards .....	7
Table 2:	Dimensions for the dimensional sketch of the movement space .....	12
Table 3:	Storage conditions.....	18
Table 4:	Operating conditions .....	18
Table 5:	Technical data .....	19
Table 6:	Sequence of error analysis .....	20

# 1 Foreword

Dear customer!



We would like to thank you for the trust you have placed in us and for purchasing our product. We have manufactured this medical product with great care.


**Please read the operating instructions carefully before using the device for the first time and always keep them within easy reach.**

These instructions for use cannot cover all possible uses of the device. For further information or in the event of problems that are not described or not described in sufficient detail in these instructions for use, please contact your specialist dealer or medical supply retailer.

# 2 General information

## 2.1 Symbols used

	<p>This warning sign indicates all instructions that are important for safety. Non-compliance can lead to accidents or injuries.</p>
	<p>Manufacturer - Shows the manufacturer of the medical device in accordance with EU Directive 2017/745 The symbol, together with the name and address of the manufacturer (i.e. the person placing the medical device on the market), must be displayed in the immediate vicinity of the symbol. appear</p>
	<p>Conformity symbol according to 2017/745 of the Medical Devices Directive</p>
	<p>Medical device - Indicates the medical device provided by the manufacturer in accordance with EU Directive 2017/745</p>
	<p>Device type B according to IEC 601-1 (Special protection against electric shock)</p>
	<p>Protection class II device, protective insulation</p>

	<p>Dispose of electrical components in accordance with legal regulations. Do not dispose of in household waste!</p>
	<p>Date of manufacture - shows the date on which the medical device was manufactured.</p>
	<p>Article number - displays the manufacturer's article number so that the medical device can be identified.</p>
	<p>Serial number - displays the manufacturer's serial number so that a specific medical device can be identified.</p>
	<p>Sales partner - indicates the company that sells the medical device locally.</p>
	<p>Temperature limit - indicates the temperature limits to which the medical device can be safely exposed.</p>
	<p>Humidity, limit - indicates the humidity range to which the medical device can be safely exposed.</p>
	<p>Air pressure, limitation - indicates the range of air pressure to which the medical device can be safely exposed.</p>
	<p>Instructions for use      Observe      or          electronic instructions for use - indicates to the user that it is necessary to follow the instructions for use.</p>

	<p>Unique identifier of a medical device - displays a carrier that contains information on a unique identifier of a medical device.</p>
	<p>Safe workload</p>
	<p>Permissible patient weight</p>
	<p>The product must be disposed of separately in the European Union. Disposal with normal household waste is not permitted.</p>
	<p>Minimum body measurements/weights of the patient</p>

## 2.2 Type plate

The type plate is attached to the head end of the trolley frame. The product can be clearly identified by the type plate.

Figure 1: Exemplary type plate

Figure 1 shows an example nameplate. Please refer to the attached type plate for the exact specifications of your product

### 2.3 Proof of standards

The following national and international standards are used for the design and verification of the product, labeling and instructions for use.

Standard	Title	Issue
DIN EN 60601-2-52	Medical electrical equipment - Part 2-52: Particular requirements for the safety of medical beds	12/2010
DIN EN 60601-1-6	Specification of suitability for use	2010
DIN EN 60601-1-2	Electromagnetic compatibility	2015
DIN EN ISO 10993	Biological evaluation of medical devices - Part 1: Assessment and verification	2010
DIN EN 1041	Provision of information by the manufacturer of a medical device	2008
DIN EN ISO 14971	Medical devices - Application of risk management to medical devices	2020

*Table 1: Verification of standards*

### 3 Safety instructions

- Before using the bed, you should read these instructions for use carefully (see Medical Devices Operator Ordinance MPBetreibV §2). It contains important information for the safe and reliable use of the device. Keep the instructions for use for future reference.
- Safety, reliability and performance are guaranteed if the following instructions are observed and the device is used properly. As the operator, you must comply with the Medical Devices Operator Ordinance (MPBetreibV).
- The Rotadorm Mobile is suitable for home care (application environment 3, 4), with a maximum patient weight of 135 kg.
- Ensure that children only have access to the bed under supervision and that no children are in the danger zone under the bed during operation.
- The bed should only be assembled by authorized personnel.
- The fuse protection on the installation side must not exceed 16A. Before connecting the appliance, please ensure that the voltage and frequency of your power supply correspond to the specifications on the rating plate.
- When choosing a location, make sure the surface is level
- Provide a suitable floor covering if the bed needs to be moved frequently. Carpets, rugs and loose floor coverings can be damaged or make it difficult to move the bed.
- Connect the mains plug firmly to the mains socket. Lay the mains connection cable on the floor. Make sure that the bed (especially when moving it) does not rest on the cable with its castors. The cable must not be routed through the mechanism of the substructure! (Danger of crushing)
- Damaged mains cables can lead to life-threatening situations. They must be replaced immediately.
- Check the mains cable for damage at regular intervals (weekly).
- Ensure that the electrical specifications of the appliance match the local conditions at the installation site.
- When the handsets are not in use, make sure that they are hanging on the bed and not placed in the bed to prevent incorrect operation that could cause damage.
- If the patient is unattended, ensure that it is set at its lowest height to make it as easy as possible to get in and out.

## 4 General product description

### 4.1 Intended use

The Rotadorm Mobile is a medical aid that can help patients to participate in their environment again and minimize physical strain for both patients and carers. The Rotadorm Mobile is the successor to the Rotadorm Care. Decisive for

This new development was based on extensive market research and the following requirements:

- Easier to divide
- Less weight of the overall bed and the individual components
- Prevention of cable breaks
- Electrical components all in the upper frame
- Weight reduction

The indications and contraindications are identical to the previous Rotadorm Care model.

### 4.2 Indication

The Rotadorm Mobile is an aid whose use is indicated:

- for people who are unable to sit up in bed and sit on the edge of the bed to get up or mobilize themselves
- to enable the patient to stand up independently
- to enable independent transfer from bed to wheelchair
- to enable and support standing up with severely restricted mobility or strain on the lumbar spine and joints
- For mobilization in cases of extreme pain symptoms during passive movement
- for extensor spasticity caused by passive movement (e.g. MS)
- for people for whom the caregiver must be physically relieved in home care in order to enable and ensure long-term home care
- for insured persons who are no longer spontaneously mobile, especially patients with QS symptoms, with a greatly increased risk of developing a pressure ulcer, or for further treatment of an existing pressure ulcer
- In cases of tetraplegia or advanced muscular dystrophy with preserved residual function of the legs, especially if regular changes of position of the upper body and legs are required
- for mobilization with apallic syndrome
- for self-mobilization after a stroke
- for seated positioning, e.g. for heart failure

The indication for the Rotadorm Mobil sit-to-stand bed applies to all illnesses or disabilities in which the patient is unable to perform the movement of standing up from lying to sitting with or without assistance due to their joints, low muscle strength or pain - or in which the caregiver is physically overloaded, e.g. due to their own illness.

### 4.3 Contraindication

The following patients are not permitted to use the Rotadorm Mobile:

- For pressure sores, especially in the hip and leg area
- In case of extreme deformation and non-weight-bearing capacity of the lower extremities

- For severe cardiovascular problems
- In the case of severe dizziness which makes it impossible to stand up independently or partially independently
- For severe anxiety

The caregiver must ensure that the user is mentally capable of handling an electrically adjustable care bed. Otherwise, all electrical functions on the bed must be switched off using the key switch and the handset must also be secured against access by the patient.

#### 4.4 Equipment features

The Rotadorm Mobile has the following electrical functions:

- Electric height adjustment
- Electric headboard adjustment
- Electric knee bend adjustment
- electrical trend/antitrend function
- Electric sit-turn function to make it easier to get out of bed

The drives for the adjustment functions consist of electromechanical linear motors with maintenance-free permanent lubrication. The drives are operated using a handset connected to the control unit via a spiral cable. The bed has four single-brake castors.

The drives and the manual switch are electrically isolated from the mains voltage and are operated with a low voltage (DC 24 V).

With a swivel seat and sit-to-stand bed, the patient's personal activity space is maintained over a long period of time. The user can move and care for themselves freely.

Depending on the equipment, the stand-up or swivel function of the lying surface can be adjusted by the service staff to suit the patient's needs. It is possible to select the direction of rotation to the left or right.

## 5 Assembly information

### 5.1 Basic information on installation

The bed should only be assembled by authorized personnel. The fuse protection on the installation side must not exceed 16A.

Before connecting the appliance, please ensure that the voltage and frequency of your power supply correspond to the specifications on the rating plate.

When choosing a location for the bed, make sure it is on a level surface. Provide a suitable floor covering if the bed needs to be moved frequently. Carpets, rugs and loose floor coverings can be damaged or make it difficult to move the bed.

Connect the mains plug firmly to the mains socket. Lay the mains connection cable on the floor. Make sure that the bed (especially when moving it) does not rest on the cable with its castors. The cable must not be routed through the mechanics of the bed base! (Danger of crushing)

**Damage to the electrical mains cable caused by driving over or clamping can have fatal consequences.**

**Before moving the bed or dismantling it for transportation, the mains connection cable must be wound up and attached to the device provided on the chassis.**

### 5.2 Dimensional sketch of the movement space

The following movement space outside the care bed must be kept free for the turning function of the lying surface. (see Figure 2 & Table 2)


Figure 2: Dimensional sketch of the movement area of the Rotadorm Mobile

Type	SI-013-1
Dimension A*)	200 cm
Dimension B*)	90 cm
Dimension C	215 cm
Dimension D	108 cm
Dimension E (head end)	40 cm
Dimension F (foot end)	50 cm

Table 2: Dimensions for the dimensional sketch of the movement space

**No objects, pieces of furniture or walls may interfere with the turning function within the specified movement range. (Danger of crushing)**

### 5.3 Dismantling the care bed

If necessary, e.g. transportation, the care beds can be dismantled with little effort as described but in reverse order. Reassembly after disassembly of the bed should be carried out by authorized personnel.

## 6 Operation

### 6.1 Operating the grab rails


Figure 3: Removable patient handle

The patient handle is used to hold or support the patient when standing up. It is only intended for use when standing up. The handle must be screwed tight using the grub screws and have latched pull catches. (see Figure 3)

**Retaining brackets must always be fitted correctly when in use.**

### 6.2 Shut-off function

All stand-up beds are equipped with a safety switch-off system - a so-called lock-out box. The integrated locking knob allows the user to lock or unlock the electrical functions directly on the box. The lock box is located underneath the lying surface.

Figure 4: Shut-off box

### 6.3 Operating the functions

All electrical functions are controlled and adjusted using the hand control. Each row of buttons is shown according to its function. (see figure 5). The first five rows of buttons are used to move the following lying surface segments in the order from top to bottom: head section, upper leg, lower leg and height. The last two rows of buttons are described in detail in the following chapters.

If necessary, an authorized specialist dealer can adjust all motor movement functions to your wishes and/or care requirements, set the individual positions of the lying surface and save them in the program sequence in the electronic control unit.


Figure 5: Patient hand control

### 6.4 Operating the Trendelenburg function

As the siphon height changes when using the Trendelenburg function (5th row of buttons), ensure that there is sufficient distance between the feet.

**Pay attention to your feet when using the Trendelenburg. Keep a sufficient distance from the bed.**

**When using the Trendelenburg function, one side must be braked, while one side is unbraked (see chapter Operating the brake rollers).**

The Trendelenburg function has an automatic locking function ("Auto-Lock"). The electric Trendelenburg function of the bed is always locked when the bed is not in use, which means that it cannot be operated using the handset. The electric Trendelenburg function can be activated using a button combination on the switch. First press the "Up" button (1) and

then press the "Down" button (2). A short beep indicates that the table has been unlocked and is now ready for use. After unlocking, the table is activated for 5 seconds. After this period, the table is automatically locked again. Locking is also signaled by a short beep.

## 6.5 Operating the rotary function

With the swivel function, the mattress base automatically moves into a suitable position for the user to get out of bed easily. (see figure 6)

Figure 6: Bed in the standing position

If the "**Stand up out**" button on the hand control is held down, all functions are run through one after the other as described below.

When the "**Stand up in**" button is pressed, the bed moves back from the stand-up position to the lying position; the head section remains slightly raised and does not move fully into the horizontal position.

Press the "**Stand up out**" button:

- The height motor moves to a predetermined level as the **starting position for turning**
- The knee bend and head section motors move simultaneously to a **predefined sitting position**, with the legs over the edge of the bed and the backrest tilted backwards by approx. 45°.
- Only then does another drive rotate the lying surface by 90° to the chassis
- After reaching the 90° position, the height of the lying surface is first corrected upwards
- then the knee bend and head section simultaneously move the lying surface into the sitting position
- The upright sitting position creates a horizontal support surface for the thigh rest, so that sitting upright is possible.
- If the "Out" button is pressed again, the bed moves to a stand-up position. This lowers the thigh support and the bed moves upwards. This makes it easier to stand up.

Press the "**Stand up in**" button:

- The height motor moves up to the preset **starting position for turning back**
- Knee bend and head section motors move simultaneously to a **predefined seat position** (see above)
- Once in the sitting position, the lying surface turns back parallel to the chassis

- The knee bend and head section motors now move into the **lying position** at the same time, with the head section remaining slightly raised and not moving back completely into the horizontal position.

**Press and hold the manual button until the turning process is complete.** (No more motor noise can be heard)

## 6.6 Operating the brake rollers

The bed has two swivel castors with double brakes. The bed should always be braked at the installation site using the castor lock. (see Figure 7)

Wheel locked **Rad gebremst**

**Rad ungebremst** Wheel unbraked

*Figure 7: Castors in braked and unbraked condition*

## 6.7 Patient lifting pole with grab handle

Our care beds are fitted with a holder for the patient trapeze bar on the left and right of the mattress base at the head end. When fitting the trapeze bar, ensure that the trapeze bar tube with the metal bolt is inserted far enough into the mounting socket so that the metal bolt is completely in the recessed groove on the mounting socket. This fixes the patient lifting pole in its position and prevents it from being swiveled beyond the lying surface. (see Figure 8)

The supplied grab handle helps the user to stand upright and can be individually adjusted to the correct height using the strap.

*Figure 8: Gallows mount at the head end*

**The trapeze bar is not suitable for therapeutic purposes.**

**The maximum load capacity of the trapeze bar is 75 kg.**

**The metal bolt on the trapeze bar must always be located in the recess groove.**

**Check the grab rail and the webbing at regular intervals for damage. Damaged parts must be replaced immediately!**

The supplied grab handle is used by the user to stand upright and can be individually adjusted to the correct height using the strap and the adjustment buckle (see Figure 9). The adjustment range is from 670mm to 870mm. (measured without mattress)

*Figure 9: Adjustable handle*

## 6.8 Operating instructions

- After assembly and before the care bed is used by a patient, check that all connections and the entire bed are secure.
- Make sure that the drives are working properly.
- A care bed that is not fully functional must be withdrawn from use immediately.
- Make sure that there are no objects, such as a wastepaper basket, side table, chair, etc., in the bed's movement area.
- To avoid the risk of injury, no parts of the patient's body should protrude above the mattress base or feet should rest on the bed frame while the adjustment functions are being operated.
- Before moving the bed, the mains plug must be disconnected from the socket to prevent damage to the electrical system.
- The bed may only be moved with patients over thresholds with a maximum height of 2 cm.
- Ensure that the duty cycle is observed. Therefore, never make long and unnecessary electrical adjustments. If the thermal fuse in the control unit is triggered once after 6 min/h, the control unit must be replaced by an authorized specialist!

**The installation of additional devices such as insulin pumps, ventilators etc. is prohibited if no potential equalization has been created beforehand.**

**Any cables from additional equipment must not be routed through the mechanics of the bed base! (Danger of crushing)**

## 6.9 Mattresses approved for use

A two-part, flame-retardant mattress in accordance with DIN 13014 and DIN 597 with a minimum density of 35 kg/m<sup>3</sup> (RG35), a compression hardness of at least 4.2 kPa, a maximum height of 12 cm, a minimum width of 88 cm and a minimum length of 197 cm (mattress and foot block together) is provided for this bed.

**The dimension from the top edge of the mattress (unloaded) to the top edge of the top wooden side rail (side rail in top position) of 22 cm must be observed for safety reasons.**

**Mattresses with a high density are only permitted if the weight of the mattress and the weight of the patient together do not exceed the safe working load of the bed.**

## 7 Ambient conditions

In accordance with DIN EN 60601-2-52, the medical device can be used in the following application environment:

### **Application environment 3:**

Long-term care in a medical setting where medical supervision is required and monitoring is provided as necessary; an ME device may be provided for a medical procedure to maintain, improve or support the patient's condition.

### **Application environment 4:**

ME device to alleviate or compensate for an injury, disability or illness in home care.

A maximum noise level of 49 dB (A) is generated during the adjustment of the electric drives.

### 7.1 Storage conditions



Storage temperature	min. -10 °C max. +50 °C	
Relative humidity	min. 20 % max. 80 %	
Air pressure (at altitude ≤ 3000 m)	min. 700 hPa max. 1060 hPa	

Table 3: Storage conditions

### 7.2 Operating conditions



Operating temperature	min. +5 °C max. +40 °C	
Relative humidity	min. 20 % max. 80 %	
Air pressure (at altitude ≤ 3000 m)	min. 700 hPa max. 1060 hPa	

Table 4: Operating conditions

## 8 Technical data Subject to change without notice!

<b>Designation</b>	<b>Type SI-013-1</b>
Nominal voltage	<b>~230 V/50Hz</b>
Rated power	<b>690 VA</b>
Device type B according to IEC 60601-1	
Protection classes	
Sound power level	<b>63 dB(A)</b>
<b>IP protection class for drive components:</b>	
Control unit	<b>IPX4</b>
Hand switch	<b>IPX4</b>
Drives	<b>IPX4</b>
Duty cycle ED 10%	<b>maximum 6 minutes/hour</b>
Safe workload	<b>175 kg</b>
Max. patient weight	<b>135 kg</b>

Dimensions of the lying surface	<b>200/210 cm x 90 cm</b>
Minimum lower leg length of the patient	<b>38 cm</b>
Bed dimensions (measured for 200 cm version)	<b>215 x 108 cm</b>
Dimensions on transport set	<b>103 x 77 x 165 cm (LxWxH)</b>
<b>Dimensions of the care bed:</b>	
Total mass incl. patient lifter and Wooden surround	<b>156.5 kg</b>
Height adjustment (measured without mattress)	<b>from 31cm to 81cm</b>
Thigh rest adjustment angle	<b>0° to 42°</b>
Adjustment angle headboard	<b>0° to 80° (with lowered leg section)</b>

Table 5: Technical data

**Repairs may only be carried out by ISKO specialist personnel or by persons authorized and trained by ISKO with comprehensive product knowledge. Any warranty and liability claims will be rejected if this provision is not complied with.**

## 9 Materials used

The care bed is manufactured as a welded tubular steel construction. The surfaces are powder-coated or galvanized. All wooden parts are either laminated or lacquered. The surfaces of this product are harmless to the skin from a health point of view.

## 10 Service and care

All household cleaners without ammonia and abrasive cleaners are permitted for cleaning the tube parts, the lying surface and the wooden parts with a damp cloth. Solvents (e.g. nitro) destroy the coating of the tubes and the lying surface!

Mechanical cleaning (e.g. scraping, sanding) or blast cleaning of the bed is not permitted. All pivot points of the moving parts, including the bearing eyes on the adjuster, are fitted with maintenance-free plain bearings and must not be oiled or greased.

## 11 Useful life of the product

With an expected average level of use in home care, the service life of the bed is 10 years. Lack of maintenance and excessive use of the product can significantly reduce the service life of the bed. The expected service life in professional nursing home use is 7 years.

## 12 Disinfection

- To ensure that the bed functions perfectly, every ISKO bed should be cleaned, disinfected and checked after each use so that it can be used again immediately.
- Improper cleaning/disinfection of the bed can be dangerous.
- Depending on the degree of soiling, we recommend cleaning the bed with a damp cloth or similar.
- Disinfectants can be used for wipe and spray disinfection in their intended concentration. (see manufacturer's instructions)

- The dilution ratio recommended by the manufacturer in the respective instructions for use must be used.

**Solvents are not permitted.**

**Abrasives or scouring pads must not be used.**

### 12.1 Specifications of the cleaning agents and disinfectants:

- The working solutions should normally be used freshly prepared.
- Do not exceed or fall below the specified concentrations.
- They must not contain any corrosive or caustic components.
- They must not contain any substances that change the surface structure or the adhesion properties of the materials.
- Lubricants must not be attacked by the cleaning agent and disinfectant.

**Under no circumstances should soap or detergents be added to the disinfectant. There is a risk of explosion and fire if alcohol-based products are used over large areas.**

**If unsuitable washing and disinfecting agents are used, damage to the surface coating may occur, for which ISKO KOCH GmbH cannot be held liable.**

## 13 Malfunctions and their rectification

	<b>Error</b>	<b>Measure</b>
1	None of the motors react to the hand switch actuations.	Check the plug connection of the mains connection cable. Check the key switch (ON/OFF)! Check the plug connection between the hand control and motor cables with the control unit.
2	One motor does not respond to the switch actuation	Check plug connection between motor cable and control unit.
3	Error number 1 or 2	Reset/initialization Press and hold the second row of buttons simultaneously (really simultaneously) until the interrupted signal tone changes to a continuous tone (after approx. 5 seconds).  Immediately after the reset (simultaneous pressing of the 2nd row of buttons), the first row of buttons is held down simultaneously for initialization until a long signal tone sounds. During this process, the the motors reach their end position.

*Table 6: Sequence of error analysis*

**Modifications, readjustments and repairs to the bed that cannot be carried out in accordance with the above instructions may only be carried out by the manufacturer directly or by a workshop authorized by the manufacturer.**

## 14 Recommended accessories

Article description	Order number
Patient handles 790 mm, coated, rigid, long	SI-031-0-25
Mattress 90 x 161 cm, incl. foot block 90 x 36 cm (200 cm)	SI-503-0
Mattress 90 x 161 cm, incl. foot block 90 x 46 cm (210 cm)	SI-503-1

## 15 Maintenance

### 15.1 Legal basis

The Medical Device Regulation (EU) 2017/745 (MDR) and national laws and regulations require operators of medical devices to ensure that the medical device is in a safe operating condition throughout its entire service life.

### 15.2 Maintenance intervals

As a requirement of the Medical Devices Operator Ordinance §4 (Maintenance), a thorough visual inspection (1), a functional test (2) and a functional test (3) must be carried out after the medical device has been in operation for two years at the latest.

(2) and a current leakage test (3) in accordance with DIN EN 62353:2015-10.

(1) Special attention must be paid to the following points during the visual inspection:

- Tight fit of all screw connections
- Mobility of the pivot points
- Check the mains supply cable for crushing or shearing points
- Checking the strain relief of the mains supply cable

(2) Special attention must be paid to the following points during the functional test:

- Function of all electrically operated movements
- Fully extend and retract all motors on the care bed (without mattress; without patient) until they switch off automatically. (Limit switches in the motors must switch off with an audible click)
- Functionality of the brakes
- Movability and function of the side rails
- Mobility of the triggers
- Checking the manual switch

**Functional tests and current leakage tests may only be carried out by ISKO specialist personnel or by persons authorized and trained by ISKO with comprehensive product knowledge.**

### 15.3 Spare parts

All spare parts for this medical device are available from ISKO KOCH GmbH, quoting the serial number, order number and article number (these can be found on the type plate attached to the medical device).

To ensure that functional reliability and any warranty claims remain valid, only original ISKO KOCH GmbH parts should be used for spare parts.

ISKO KOCH GmbH

Egerländer Street 28

95448 Bayreuth

Phone: +49(0)921/150845-0 (Monday - Thursday 8:00 - 17:00 & Friday 8:00 - 16:00) Fax:

+49(0)921/150845-45

E-mail: [info@isko-koch.de](mailto:info@isko-koch.de)

## 15.4 Notes on documentation

According to the Medical Device Operator Ordinance and the Medical Device Regulation (EU) 2017/745 (MDR), there is a documentation obligation for

- Maintenance
- Incidents / near misses

If unusual hazards for the product are foreseeable at the installation site of the care bed (supply line is on the floor; children playing; pets; ...), the electrical cables in particular must be constantly checked and suitable measures taken to avert hazards.

## 16 Re-use

Before each reuse of the care bed, a thorough visual and functional check of all electrically operated functions as well as a current leakage test in accordance with DIN EN 62353:2015-10 must be carried out as described in the **maintenance intervals** section. The points on **service and care & maintenance** specified in the operating instructions must be observed when cleaning the bed.

## 17 Waste disposal

### 17.1 Disposal of the device

The appliance and accessories, if present, should be disposed of in an environmentally friendly manner and in accordance with legal regulations. Please comply with the applicable **waste separation regulations!** If you are still unsure about this, please contact your local authority or waste disposal company.

### 17.2 Disposal of electrical components

\*if electrical components are included in the medical device

In accordance with Directive 2012/19/EU - WEEE2, this medical device is classified as electrical equipment. All electrical components are free of non-authorized substances classified as harmful in accordance with RoHS II Directive 2011/65/EU. In addition, replaced electrical components must be disposed of in accordance with the European directives (see Directive 2012/19/EU - WEEE2).

### 17.3 Disposal of the packaging

The EU Waste Framework Directive 2008/98/EC is decisive for the handling of packaging disposal. Reusable materials must be recycled in accordance with national regulations.

## 18 EC Declaration of Conformity

As the manufacturer, we declare under our sole responsibility that our sit-to-stand beds comply with the essential requirements of the EC Directive for

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