Stehgeräte

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Instruction for use

Stand: 09/2022  
(Rev. 4.0)

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**Version: 09/2022  
(Rev. 4.0)**

## Table of contents

1	Foreword .....	6
2	General information .....	7
2.1	Symbols used .....	7
2.2	Type plate .....	7
2.3	Proof of standards .....	10
3	Safety instructions .....	11
4	General product description .....	12
4.1	Intended use .....	12
4.2	Indication .....	12
4.3	Contraindication .....	12
4.4	Equipment features .....	12
4.5	Scope of delivery .....	12
5	Assembly information .....	13
5.1	Basic information on installation .....	13
6	Service .....	13
6.1	Body size .....	13
6.2	Back pad or hip pad .....	14
6.3	Foot and heel support .....	15
6.4	Knee pads .....	16
6.5	Upper body pad .....	16
6.5.1	Adjusting the backrest cushion .....	16
6.5.2	Lateral armrest .....	16
6.5.3	Adjusting the width of the lateral hip supports. ....	16
6.5.4	Attaching and adjusting the headrest .....	16
6.6	Breast pad .....	17
6.7	General instructions for use when standing up. ....	18
6.8	Operating instructions .....	20
6.9	Behavior in the event of a power failure .....	20
6.10	Operating the brake rollers .....	21
7	Ambient conditions .....	22
7.1	Storage conditions .....	22
7.2	Operating conditions .....	22
8	Technical data subject to change! .....	23
9	Materials used .....	23
10	Service and care .....	23

11	Useful life of the product .....	23
12	Disinfection.....	24
12.1	Specifications of cleaning agents and disinfectants:.....	24
13	Malfunctions and their rectification.....	24
14	Recommended accessories .....	25
15	Maintenance. ....	26
15.1	Legal basis. ....	26
15.2	Maintenance intervals.....	26
15.3	Spare parts .....	27
15.4	Notes on documentation .....	27
16	Re-use.....	27
17	Waste disposal.....	27
17.1	Disposal of the device .....	27
17.2	Disposal of electrical components.....	27
17.3	Disposal of packaging .....	27
18	EC Declaration of Conformity .....	28

## List of illustrations

Figure 1: Exemplary type plate .....	10
Figure 2: Adjusting the table height.....	15
Figure 3: Clamping lever for table height adjustment .....	15
Figure 4: Position of the hip pad in the buttocks.....	16
Figure 5: Adjusting the heel switches and toe fixation.....	16
Figure 6: Adjusting the knee pads.....	17
Figure 7: Upper body pad.....	18
Figure 8: Breast pad .....	18
Figure 9: Position of the patient restraint belt in the buttocks when the patient is in the wheelchair	19
Figure 10: Hand control .....	20
Figure 11: Electric belt retractor .....	20
Figure 12: Example castors in braked and unbraked condition .....	22

## List of tables

Table 1: Symbols used.....	10
Table 2: Verification of standards .....	11
Table 3: Storage conditions.....	23
Table 4: Operating conditions .....	23
Table 5: Technical data .....	24
Table 6: Malfunctions and their rectification .....	25
Table 7: Recommended accessories .....	26

## 1 Foreword

Dear customer!



We would like to thank you for the trust you have placed in us and for purchasing our product. We have manufactured this medical product with great care.


**Please read the operating instructions carefully before using the device for the first time and always keep them within easy reach.**

These instructions for use cannot cover all possible uses of the device. For further information or in the event of problems that are not described or not described in sufficient detail in these instructions for use, please contact your specialist dealer or medical supply retailer.

## 2 General information

### 2.1 Symbols used

	<p>This warning sign indicates all instructions that are important for safety. Non-compliance can lead to accidents or injuries.</p>
	<p>Manufacturer - Shows the manufacturer of the medical device in accordance with EU Directive 2017/745 The symbol, together with the name and address of the manufacturer (i.e. the person placing the medical device on the market), must be placed in the immediate vicinity of the symbol. appear</p>
	<p>Conformity symbol according to 2017/745 of the Medical Devices Directive</p>
	<p>Medical device - Indicates the medical device provided by the manufacturer in accordance with EU Directive 2017/745</p>
	<p>Device type B according to IEC 601-1 (Special protection against electric shock)</p>
	<p>Protection class II device, protective insulation</p>
	<p>Dispose of electrical components in accordance with legal regulations. Do not dispose of in household waste!</p>
	<p>Date of manufacture - shows the date on which the medical device was manufactured.</p>
	<p>Article number - displays the manufacturer's article number so that the medical device can be identified.</p>

	<p>Serial number - displays the manufacturer's serial number so that a specific medical device can be identified.</p>
	<p>Sales partner - indicates the company that sells the medical device locally.</p>
	<p>Temperature limit - indicates the temperature limits to which the medical device can be safely exposed.</p>
	<p>Humidity, limit - indicates the humidity range to which the medical device can be safely exposed.</p>
	<p>Air pressure, limitation - indicates the range of air pressure to which the medical device can be safely exposed.</p>
	<p>Instructions for use Observe or electronic instructions for use - indicates to the user that it is necessary to follow the instructions for use.</p>
	<p>Unique identifier of a medical device - displays a carrier that contains information on a unique identifier of a medical device.</p>
	<p>Safe workload</p>
	<p>Permissible patient weight</p>


	<p>The product must be disposed of separately in the European Union. Disposal with normal household waste is not permitted.</p>
	<p>Minimum body measurements/weights of the patient</p>

Table 1: Symbols used

## 2.2 Type plate

The type plate is attached to the head end of the trolley frame. The product can be clearly identified by the type plate.


Figure 1: Exemplary type plate

Figure 1 shows an example nameplate. Please refer to the attached type plate for the exact specifications of your product

### 2.3 Proof of standards

The following national and international standards are used for the design and verification of the product, labeling and instructions for use.

<b>Standard</b>	<b>Title</b>	<b>Issue</b>
DIN EN 60601-1-6	Specification of suitability for use	2010
DIN EN 60601-1-2	Electromagnetic compatibility	2015
DIN EN ISO 10993	Biological evaluation of medical devices - Part 1: Assessment and verification	2010
DIN EN 1041	Provision of information by the manufacturer of a medical device	2008
DIN EN ISO 14971	Medical devices - Application of risk management to medical devices	2020

*Table 2: Verification of standards*

### 3 Safety instructions

- Before using the standing device, you should read these instructions for use carefully (see Medical Devices Operator Ordinance MPBetreibV §2). It contains important information for the safe and reliable use of the device. Keep the instructions for use in a safe place for future reference.
- Safety, reliability and performance are guaranteed if the following instructions are observed and the device is used properly. As the operator, you must comply with the Medical Devices Operator Ordinance (MPBetreibV).
- The device may only be operated by medically trained and instructed personnel at the start of upright training. Operators must be trained in the use of the device and their training and level of knowledge must enable them to assess the risks for the patient. The trained personnel are responsible for the assessment for subsequent independent use by the user (patient), whereby an assistant must remain within call or sight.
- All clamp and screw connections must be checked regularly
- During use, 1 or 2 assistants must be available to support the user.
- Only use aids on a firm, level and horizontal surface
- Do not use the aid outdoors. The device may only be used in protected domestic and clinical areas.
- Leaning out far poses a risk of falling
- **The aid must not be used as a means of transportation**
- Only park the aid with the wheels locked
- In the event of a fault, the standing appliance must be withdrawn from use immediately
- Parts of the standing appliance can catch fire if handled improperly
- When standing, the user must always wear sturdy shoes; barefoot use is not permitted. Unless expressly instructed by the therapist
- Ensure that children only have access to the standing appliance under supervision.
- The standing device should only be set up by authorized personnel.
- The fuse protection on the installation side must not exceed 16A. Before connecting the appliance, please ensure that the voltage and frequency of your power supply correspond to the specifications on the rating plate.
- When choosing a location, make sure the surface is level
- Provide a suitable floor covering if the standing appliance has to be moved frequently. Carpets, rugs and loose floor coverings can be damaged or make it difficult to move the appliance.
- Connect the mains plug firmly to the mains socket. Lay the mains connection cable on the floor. Make sure that the upright appliance (especially when moving it) does not stand on the cable with its castors. The cable must not be routed through the mechanics of the substructure! (Danger of crushing)
- Damaged mains cables can lead to life-threatening situations. They must be replaced immediately.
- Check the mains cable for damage at regular intervals (weekly).
- Ensure that the electrical specifications of the appliance match the local conditions at the installation site.
- When the handsets are not in use, ensure that they are attached to the standing appliance to prevent incorrect operation that could cause damage.

## 4 General product description

### 4.1 Intended use

The Levacare standing device is used exclusively for standing training for people who are unable to walk or have difficulty walking for independent use with an assistant. Standing training is used in inpatient and outpatient rehabilitation as well as in the home environment. The design of our standing devices

enables independent, stable and safe standing. An electric belt retractor brings the patient out of the wheelchair and into a secure standing position.

### 4.2 Indication

People between 140 and 200 cm tall with various physical and/or psychological limitations can stand safely, firmly and without fatigue every day. Daily standing stabilizes the circulation, prevents bone decalcification (osteoporosis), corrects contractures of the hip and knee joints, stimulates bowel and bladder functions as well as breathing and kidney activity. Targeted promotion of head control is optimized by the precise adjustment options at all levels.

### 4.3 Contraindication

The following patients are not permitted to use the standing device:

- For pressure sores, especially in the hip and leg area
- In case of extreme deformation and non-weight-bearing capacity of the lower extremities
- For severe cardiovascular problems
- In the case of severe dizziness which makes it impossible to stand up independently or partially independently
- For severe anxiety

### 4.4 Equipment features

The Levacare standing device has the following features:

- Table cut-out with ergonomic belly pad
- Foot plate with heel cups and toe fixation
- Double swivel castors with total locking
- Table height adjustment with gas spring support
- Knee pads adjustable in height and depth, as well as with horizontal and angle adjustment
- Stainless steel handrails height-adjustable in parallel
- Electric belt retractor system available
- Simple maintenance and repair of the electrical or mechanical belt retractor system

### 4.5 Scope of delivery

- Pre-installed stable base frame with bar rails adjustable from 84 to 107 cm
- Support table with height adjustment from 100 to 130 cm
- Height and depth-adjustable knee pads
- Heel cups with Velcro fasteners

#### **Optional:**

- Stand-up aid with electric automatic belt tensioner

## 5 Assembly information

### 5.1 Basic information on installation

The standing appliance should only be set up by authorized personnel. The fuse protection on the installation side must not exceed 16A.

Before connecting the appliance, please ensure that the voltage and frequency of your power supply correspond to the specifications on the rating plate.

Ensure that the floor is level when selecting a location for the standing device. Provide a suitable floor covering if the standing device has to be moved frequently. Carpets, rugs and loose floor coverings can be damaged or make it difficult to move.

Connect the mains plug firmly to the mains socket. Lay the mains connection cable on the floor. Make sure that the upright appliance (especially when moving it) does not rest on the cable with its castors.

**Damage to the electrical mains cable caused by driving over or clamping can have fatal consequences.**

## 6 Operation

All presettings may only be made on the empty standing device without the patient standing in it. The settings must then be checked again with the patient standing in the aid and corrected if necessary. An assistant is required for the adjustment work. Due to the large number of setting options for the standing device, some of the settings described may not be possible for some users or under certain circumstances.

**Check that all connections are tight after adjustment work!**

**The cushion of the table top must never be set below the top edge of the pool!**

### 6.1 Body size

Adjustment to the patient's height is made via the extension height of the table tubes. The adjustment is correct when the therapy table is at arm height or slightly higher if required for a standing patient (see Figure 2). The table height adjustment is supported by two built-in gas springs.

Loosen the clamping levers on the standpipe and pull it out as required (see Figure 3).

**The clamping lever must always be tightened when using the standing device!**


Figure 2: Adjusting the table height

Figure 3: Clamping lever for table height adjustment

## 6.2 Back pad or hip pad

**An assistant is required to use the back cushions!**

The height of the back cushion or hip support is adjusted by loosening the clamping levers on the main tubes and the depth is adjusted by loosening the clamping levers on the rear cross tube. Adjustment to the width is made using the catches on the rear cross tube. The pads are positioned above the buttocks in the area of the lumbar vertebrae when the patient is standing (see Figure 4). They should have a stabilizing effect in conjunction with the heel support and knee pad.

**The clamping lever must always be tightened when using the standing device!**

Figure 4: Position of the hip pad in the buttocks

### 6.3 Foot and heel holder

Adjust the heel holders by repositioning them on the perforated grid plate so that the user has a firm hold (4 positions are possible).

The patient can get in by opening the Velcro straps on the harness. The straps are closed tightly over the instep. The Velcro straps serve as toe fixation. (see Figure 5)

**Check that the foot is firmly attached!**

Figure 5: Adjusting the heel switches and toe fixation

## 6.4 Knee pads

The height and depth of the knee pads are adjusted by loosening the vertical and horizontal tube catches. The knee pads are correctly adjusted when they are positioned below the knee on the shin. The knee pads should first be preset and corrected if necessary when the patient is standing. They stabilize the patient's stance in conjunction with the heel support and back pad. (see Figure 6)

*Figure 6: Adjusting the knee pads*

## 6.5 Upper body pad

(Optional as accessory item number KR-630-X)

### 6.5.1 Adjusting the backrest cushion

The back pad supports the middle of the back at shoulder blade level. To attach the pad, turn the knob on the holder a full turn and pull it out. Slide the pad holder into the bracket with the holes facing outwards. Lower the tube into the holder and release the knob as soon as the breast pad is at the desired height. Release the knob so that the indexing pin integrated into the knob engages in the nearest hole in the tube. Tighten the knob to fix the back pad in position.

### 6.5.2 Lateral armrest

The chest pad is supplemented by side pads. Two width-adjustable pads can be attached to the chest pad, which are positioned on both sides of the chest at the height of the armpits. The height of the pads depends on the height adjustment of the chest pad and should be adjusted so that there is a distance of three fingers between the top of the side chest pads and the armpit. See previous instructions for adjusting the height of the breast pads.

### 6.5.3 Adjusting the width of the lateral hip supports

The hip supports can be individually adjusted in width to center the thorax. Their position is held in place with locking pins. Pull out the plungers and slide the arm supports along the horizontal tube to the desired position. Snap the plungers into the nearest holes and make sure that they are locked in the desired position. Ensure that the supports are positioned symmetrically.

### 6.5.4 Attaching and adjusting the headrest

The headrest supports the head at the level of the back of the head. To attach the headrest, pull out the locking pin on the chest support. Slide the tube of the headrest with the holes facing

into the holder on the outside. When you have reached the desired height, release the indexing plunger so that it snaps into the nearest hole in the headrest tube.


The depth of the headrest should be adjusted so that the head is supported but not pushed forward. To adjust the depth of the headrest, pull the plunger on the top of the headrest upwards and slide the headrest in or out to the desired position. When you have reached the desired position, release the plunger and make sure that it engages in one of the holes in the headrest.

*Figure 7: Upper body pad*

## 6.6 Breast pad

The chest support supports the patient when standing. After centering and aligning the chest support on the table, it can be attached to the table of the standing device using two screw clamps. (see Figure 8)

*Figure 8: Breast pad*

## 6.7 General instructions for use when standing up

**Put on the patient harness before approaching the standing device.**

The patient belt should be positioned so that its lower edge is clamped between the buttocks and the wheelchair seat. Then fasten the abdominal belt (see Figure 9).

Swivel the legrests or footrest to the side. Place the feet in the foot guides. The heel sections can be removed for this purpose. Drive the wheelchair so close that the front wheels touch the recesses provided in the footplate. Secure the wheelchair against rolling away. (see Figure 5)

*Figure 9: Position of the patient restraint belt in the buttocks when the patient is in the wheelchair*

Once the electric belt retractor (see Figure 11) has been connected to the mains using the mains plug, the device is ready for operation. The retractor can be operated with the hand switch (see Figure 10) or optionally with the suction/blow hose. The drive is automatically switched off in the two end bearings of the belts. This prevents the patient from becoming trapped under any circumstances.

The additional forced switch-off is triggered as long as the switch-off button (pos. 3 in Figure 11) is pressed. The device can only be operated again once it has been released. Each upward or downward movement of the straps is effected by pressing the respective button on the hand control. The buttons are marked with arrows according to their function. If suction is applied to the mouthpiece of the hose, the straps are moved up and down when the mouthpiece is blown. The air hose must always be laid and stored in such a way that it cannot be accidentally operated by the user, an assistant or a third party. Malfunctions can occur if the air hose is kinked or crushed. If the air switch is not required, the air hose can be disconnected from the connector.

Figure 10: Manual switch

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

Figure 11: Electric belt retractor

## 6.8 Operating instructions

Individual adjustment is required before using the standing device with electric belt retractor. It is particularly important for the safe functioning of all equipment that the table height and the switch-off button for the forced shutdown are set correctly. The table height must be adjusted by loosening the clamping levers so that the cushion of the table top is between the sternum and the upper edge of the pelvis in the standing position according to personal requirements. (see chapter 6.1)

The cushion of the table top must never be set below the top edge of the pool!

**Once the adjustment is complete, tighten the clamping levers again.**

To set the shut-off button of the forced shut-off 3, loosen the lock nut 6 and turn the shut-off button on the thread of the switch rod 7. The shut-off button is set correctly if the motor is not yet shut off up to the desired stop position, but shut-off is imminent.

**Once the setting has been made, the switch-off button must be secured again with the lock nut.**

When replacing or changing the belt retractor, make sure that no one is standing by the appliance when this is being done. The appliance must also be disconnected from the power supply.

Now remove the two locking screws 9 on the two guide rods. Approach the appliance from behind and support the belt retractor with your forearms to prevent uncontrolled sliding when pulling the locking bolts. You can now slide the belt retractor down the retaining rods and remove it.

Attach in reverse order. The appliance is correctly fitted when the locking bolts are correctly engaged and the locking screws 9 are in place.

**All functions of the belt retractor are switched off while the forced switch-off button is pressed. Further operation is only possible when the switch-off button is released; if necessary, press it back by hand. When unwinding the belts, make sure that they are tensioned to prevent the belts from becoming entangled.**

## 6.9 Behavior in the event of a power failure

If the power fails during use of the electric belt retractor, the user must call an assistant to help the patient back into the wheelchair.

The following description assumes that the assistant is right-handed. For a left-handed person, proceed in reverse.

First, the wheelchair is turned around the left front wheel just enough so that the assistant can stand behind the patient, but it is still possible to place the patient in the wheelchair. The assistant stands behind the user, opens the Velcro fastener, grasps the patient's chest with the right arm under the patient's right armpit and opens the left buckle of the patient restraint belt with the left hand (if necessary, relieve the patient restraint belt by applying a little pressure via the body buckle). The assistant's left hand now grasps under the patient's buttocks to support them as they slide back into the wheelchair.

## 6.10 Operating the brake rollers

The standing device has four swivel castors with individual brakes. The standing device should always be braked at the installation site and during use using the castor lock. (see Figure 12)


Figure 12: Exemplary rollers in braked and unbraked condition

**Standing training with unbraked castors can lead to injuries to the patient.**

## 7 Ambient conditions

### 7.1 Storage conditions




Storage temperature	min. -10 °C max. +50 °C	
Relative humidity	min. 20 % max. 80 %	
Air pressure (at altitude ≤ 3000 m)	min. 700 hPa max. 1060 hPa	

Table 3: Storage conditions

### 7.2 Operating conditions



Operating temperature	min. +5 °C max. +40 °C	
Relative humidity	min. 20 % max. 80 %	
Air pressure (at altitude ≤ 3000 m)	min. 700 hPa max. 1060 hPa	

Table 4: Operating conditions

## 8 Technical data subject to change!

Designation	KR-100-B	KR-111-B
Nominal voltage		~230 V/50Hz
Rated power		690 VA
Device type B according to IEC 601-1		
Protection class		
Sound power level		63 dB(A)
IP protection class Hand switch, control box, Motor		IP X4
IP protection class Standing device lifting mechanism		IP X1
Duty cycle ED 10%		maximum 6 Minutes/hour
max. patient weight	120 kg	120 kg
Patient size	140 - 200 cm	140 - 200 cm
Empty weight	approx. 52 kg	approx. 65 kg
Width	77 cm	77 cm
Length	118 cm	118 cm
Table height	100 - 130 cm	100 - 130 cm

Table 5: Technical data

**Repairs may only be carried out by ISKO specialist personnel or by persons authorized and trained by ISKO with comprehensive product knowledge. Any warranty and liability claims will be rejected if this provision is not complied with.**

## 9 Materials used

The upright unit is manufactured as a welded tubular steel construction. The surfaces are powder-coated or galvanized. All wooden parts are either laminated or lacquered. The surfaces of this product are harmless to the skin from a health point of view.

## 10 Service and care

All household cleaners without ammonia and abrasive cleaners are permitted for cleaning the tube parts, the lying surface and the wooden parts with a damp cloth. Solvents (e.g. nitro) destroy the coating of the tubes and the lying surface!

Mechanical cleaning (e.g. scraping, sanding) or blast cleaning of the standing device is not permitted. All pivot points of the moving parts, including the bearing eyes on the adjuster, are fitted with maintenance-free plain bearings and must not be oiled or greased.

## 11 Useful life of the product

With an expected average level of use in home care, the service life of the standing device is 10 years. Lack of maintenance and excessive strain on the product can significantly reduce the service life of the standing device. The expected service life in professional nursing home use is 7 years.

## 12 Disinfection

- To ensure that the standing device functions perfectly, every ISKO standing device should be cleaned, disinfected and checked after each use so that it can be used again immediately.
- Improper cleaning/disinfection of the standing appliance can be dangerous.
- Depending on the degree of soiling, we recommend cleaning the standing appliance with a damp cloth or similar.
- Disinfectants can be used for wipe and spray disinfection in their intended concentration. (see manufacturer's instructions)
- The dilution ratio recommended by the manufacturer in the respective instructions for use must be used.

**Solvents are not permitted.**

**Abrasives or scouring pads must not be used.**

### 12.1 Specifications of the cleaning agents and disinfectants:

- The working solutions should normally be used freshly prepared.
- Do not exceed or fall below the specified concentrations.
- They must not contain any corrosive or caustic components.
- They must not contain any substances that change the surface structure or the adhesion properties of the materials.
- Lubricants must not be attacked by the cleaning agent and disinfectant.

**Under no circumstances should soap or detergents be added to the disinfectant. There is a risk of explosion and fire if alcohol-based products are used over large areas.**

**The use of unsuitable detergents and disinfectants may cause damage to the surface coating, for which ISKO KOCH GmbH cannot be held liable.**

## 13 Malfunctions and their rectification

	<b>Error</b>	<b>Measure</b>
1	None of the motors react to the manual button actuations.	Check the plug connection of the mains connection cable. Check the key switch (ON/OFF)! Check the plug connection between the hand control and motor cables with the control unit.
2	One motor does not respond to the switch actuation	Check plug connection between motor cable and control unit.

*Table 6: Malfunctions and their rectification*

**Modifications, readjustments and repairs to the standing device that cannot be carried out in accordance with the above instructions may only be carried out by the manufacturer directly or by a workshop authorized by the manufacturer.**

## 14 Recommended accessories

Article description	Order number
Support bar with hip supports	KR-620-X
Straight upper body support with head and arm support	KR-631-X
Breast pad 15 cm, is attached to the table with screw clamps	KR-744-X
Breast pad 25 cm, is attached to the table with screw clamps	KR-745-X
Suction/blower control for electric belt retractor	KR-601-SB
Table top horizontally adjustable	KR-710-X
Foot plate horizontally adjustable	KR-720-X
Table edge	KR-602-X
Electric belt retractor	KR-601-X
Patient restraint belt, size S, standard	KR-625-X
Patient restraint belt, size M, standard	KR-624-X
Patient restraint belt, size L, standard	KR-623-X
Patient restraint belt, size S, with base	KR-608-X
Patient restraint belt, size M, with base	KR-607-X
Patient restraint belt, size L, with base	KR-606-X
Patient restraint belt, size XL, with base	KR-609-X

*Table 7: Recommended accessories*

## 15 Maintenance

### 15.1 Legal basis

The Medical Device Regulation (EU) 2017/745 (MDR) and national laws and regulations require operators of medical devices to ensure that the medical device is in a safe operating condition throughout its entire service life.

### 15.2 Maintenance intervals

As a requirement of the Medical Devices Operator Ordinance §4 (Maintenance), a thorough visual inspection (1), a functional test (2) and a functional test (3) must be carried out after the medical device has been in operation for two years at the latest.

(2) and a current leakage test (3) in accordance with DIN EN 62353:2015-10.

(1) Special attention must be paid to the following points during the visual inspection:

- Tight fit of all screw connections
- Mobility of the pivot points
- Check the mains supply cable for crushing or shearing points
- Checking the strain relief of the mains supply cable

(2) Special attention must be paid to the following points during the functional test:

- Function of all electrically operated movements
- Fully extend and retract all motors on the standing device (without patient) until they switch off automatically. (Limit switches in the motors must switch off with an audible click)
- Functionality of the brakes
- Mobility of the triggers
- Checking the manual switch

**Functional tests and current leakage tests may only be carried out by ISKO specialist personnel or by persons authorized and trained by ISKO with comprehensive product knowledge.**

### 15.3 Spare parts

All spare parts for this medical device are available from ISKO KOCH GmbH, stating the serial number, order number and article number (these can be found on the type plate attached to the medical device).

To ensure that functional reliability and any warranty claims remain valid, only original ISKO KOCH GmbH parts should be used for spare parts.

ISKO KOCH GmbH

Egerländer Street 28

95448 Bayreuth

Phone: +49(0)921/150845-0 (Monday - Thursday 8:00 - 17:00 & Friday 8:00 - 16:00) Fax:

+49(0)921/150845-45

E-mail: [info@isko-koch.de](mailto:info@isko-koch.de)

## 15.4 Notes on documentation

According to the Medical Device Operator Ordinance and Medical Device Regulation (EU) 2017/745 (MDR), there is a documentation obligation for:

- Maintenance
- Incidents / near misses

If unusual hazards for the product are foreseeable at the installation site of the standing appliance (supply line is on the ground; children playing; pets; ...), the electrical cables in particular must be constantly checked and suitable measures taken to avert hazards.

## 16 Re-use

Before using the standing device again, a thorough visual and functional check of all electrically operated functions as well as a current leakage test in accordance with DIN EN 62353:2015-10 must be carried out as described under **Maintenance intervals**. The points on **service and care & maintenance** specified in the operating instructions must be observed when cleaning the standing device.

## 17 Waste disposal

### 17.1 Disposal of the device

The appliance and accessories, if present, should be disposed of in an environmentally friendly manner and in accordance with legal regulations. Please comply with the applicable **waste separation regulations!** If you are still unsure about this, please contact your local authority or waste disposal company.

### 17.2 Disposal of electrical components

\*if electrical components are included in the medical device

In accordance with Directive 2012/19/EU - WEEE2, this medical device is classified as electrical equipment. All electrical components are free of non-authorized substances classified as harmful in accordance with RoHS II Directive 2011/65/EU. In addition, replaced electrical components must be disposed of in accordance with the European directives (see Directive 2012/19/EU - WEEE2).

### 17.3 Disposal of the packaging

The EU Waste Framework Directive 2008/98/EC is decisive for the handling of packaging disposal. Reusable materials must be recycled in accordance with national regulations.

## 18 EC Declaration of Conformity

As the manufacturer, we declare under our sole responsibility that our standing appliances comply with the essential requirements of the EC Directive for

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