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Stand: 05/2022
(Rev. 2.0)

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
**Status: 05/2022
(Rev. 2.0)**

Table of contents

1	Foreword.....	4
2	General information.....	4
2.1	Symbols used	4
2.2	Type plate.....	6
2.3	Proof of standards.....	7
3	Safety instructions.....	7
4	General product description	7
4.1	Purpose/indication	7
4.2	Contraindication.....	7
5	Conversion process	8
6	Ambient conditions	8
7	Technical data subject to change!.....	9
8	Materials used.....	9
9	Service and maintenance	9
10	Useful life of the product	9
11	Maintenance	10
11.1	Legal basis	10
11.2	Maintenance intervals.....	10
11.3	Spare parts	10
11.4	Notes on documentation	10
12	Re-use	10
13	Waste disposal	11
13.1	Disposal of the device	11
13.2	Disposal of electrical components.....	11
13.3	Disposal of the packaging.....	11
14	EC Declaration of Conformity	12

List of illustrations

Figure 1: Exemplary type plate	6
Figure 2: Transfer process from wheelchair to care bed as an example	8

List of tables

Table 1: Symbols used	6
Table 2: Verification of standards	7
Table 3: Operating conditions	8
Table 4: Technical data	9

1 Foreword

Dear customer!



We would like to thank you for the trust you have placed in us and for purchasing our product. We have manufactured this medical product with great care.


Please read the operating instructions carefully before using the device for the first time and always keep them within easy reach.

These instructions for use cannot cover all possible uses of the device. For further information or in the event of problems that are not described or not described in sufficient detail in these instructions for use, please contact your specialist dealer or medical supply retailer.

2 General information

2.1 Symbols used

	<p>This warning sign indicates all instructions that are important for safety. Non-compliance can lead to accidents or injuries.</p>
	<p>Manufacturer - Shows the manufacturer of the medical device in accordance with EU Directive 2017/745 The symbol, together with the name and address of the manufacturer (i.e. the person placing the medical device on the market), must be displayed in the immediate vicinity of the symbol. appear</p>
	<p>Conformity symbol according to 2017/745 of the Medical Devices Directive</p>
	<p>Medical device - Indicates the medical device provided by the manufacturer in accordance with EU Directive 2017/745</p>
	<p>Device type B according to IEC 601-1 (Special protection against electric shock)</p>
	<p>Protection class II device, protective insulation</p>

	<p>Dispose of electrical components in accordance with legal regulations. Do not dispose of in household waste!</p>
	<p>Date of manufacture - shows the date on which the medical device was manufactured.</p>
	<p>Article number - displays the manufacturer's article number so that the medical device can be identified.</p>
	<p>Serial number - displays the manufacturer's serial number so that a specific medical device can be identified.</p>
	<p>Sales partner - indicates the company that sells the medical device locally.</p>
	<p>Temperature limit - indicates the temperature limits to which the medical device can be safely exposed.</p>
	<p>Humidity, limit - indicates the humidity range to which the medical device can be safely exposed.</p>
	<p>Air pressure, limitation - indicates the range of air pressure to which the medical device can be safely exposed.</p>
	<p>Instructions for use Observe or electronic instructions for use - indicates to the user that it is necessary to follow the instructions for use.</p>

	Unique identifier of a medical device - displays a carrier that contains information on a unique identifier of a medical device.
	Safe workload
	Permissible patient weight
	Minimum body measurements/weights of the patient

Table 1: Symbols used

2.2 Type plate

The type plate is attached to the head end of the trolley frame. The type plate allows the product to be clearly identified.

Figure 1: Exemplary type plate

Figure 1 shows an example type plate. Please refer to the attached type plate for the exact specifications of your product.

2.3 Proof of standards

The following national and international standards are used for the design and verification of the product, labeling and instructions for use.

Standard	Title	Issue
DIN EN 60601-1-6	Specification of suitability for use	2010
DIN EN ISO 10993	Biological evaluation of medical devices - Part 1: Assessment and verification	2010
DIN EN 1041	Provision of information by the manufacturer of a medical device	2008
DIN EN ISO 14971	Medical devices - Application of risk management to medical devices	2020

Table 2: Verification of standards

3 Safety instructions

- Before using the product, you should read these instructions for use carefully (see Medical Devices Operator Ordinance MPBetreibV §2). It contains important information for the safe and reliable use of the device. Keep the instructions for use in a safe place for future reference.
- Safety, reliability and performance are guaranteed if the following instructions are observed and the device is used properly. As the operator, you must comply with the Medical Devices Operator Ordinance (MPBetreibV).

4 General product description

4.1 Purpose/indication

The sliding board enables patients/disabled persons who still have sufficiently strong arm muscles and trunk stability to change position without the need for an assistant. The slide board is a board with a very smooth surface in order to minimize frictional resistance and thus prevent slipping.

to enable the patient to glide over it easily. It is used to transfer the patient to a wheelchair or care bed, for example. The patient can change position independently or, if necessary, with assistance. In general, instruction and a check of the user's physical condition by specialist therapeutic staff are necessary for use.

4.2 Contraindication

The following patients are not permitted to use the sliding board:

- Decubitus in the buttock area
- Temporary or permanent dizziness
- Impaired cognitive perception

5 Conversion process

Figure 2: Transfer process from wheelchair to care bed as an example

1. Slide the sliding board under the buttocks
2. Then grasp the end with one hand
3. Now pull your body over the slide board

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Table 4: Technical data

A "D" is added after the order number if the thickness changes from 8 mm to 16 mm. Without this "D", the 8 mm variant is indicated by the order number.

8 Materials used

The surface of the sliding board is made of melamine resin Eurodecor. The surfaces of this product are unthinkable for the skin from a health point of view.

9 Service and care

All household cleaners without ammonia and abrasive cleaners are permitted for cleaning the sliding board with a damp cloth.

Please note: Solvents (e.g. nitro) destroy the coating, mechanical cleaning (e.g. sanding, scraping) is not permitted. Do not use sharp-edged tools (e.g. knives, metal scrapers) for cleaning.

The slide board is intended for use in domestic or protected outdoor areas.

10 Useful life of the product

With an expected average level of use in home care, the service life of the sliding board is 10 years. Lack of maintenance and excessive use of the product can significantly reduce the service life of the sliding board. The expected service life in professional care home use is 7 years.

11 Maintenance

11.1 Legal basis

The Medical Device Regulation (EU) 2017/745 (MDR) and national laws and regulations require operators of medical devices to ensure that the medical device is in a safe operating condition throughout its entire service life.

11.2 Maintenance intervals

As a requirement of the Medical Devices Operator Ordinance §4 (Maintenance), a thorough visual inspection (1) must be carried out after the medical device has been in operation for two years at the latest.

(1) Special attention must be paid to the following points during the visual inspection:

- Cracks
- Surface injuries
- Fracture points

11.3 Spare parts

All spare parts for this medical device are available from ISKO KOCH GmbH, quoting the serial number, order number and article number (these can be found on the type plate attached to the medical device).

To ensure that functional reliability and any warranty claims remain valid, only original ISKO KOCH GmbH parts should be used for spare parts.

ISKO KOCH GmbH

Egerländer Street 28

95448 Bayreuth

Phone: +49(0)921/150845-0 (Monday - Thursday 8:00 - 17:00 & Friday 8:00 - 16:00) Fax:

+49(0)921/150845-45

E-mail: info@isko-koch.de

11.4 Notes on documentation

According to the Medical Device Operator Ordinance and Medical Device Regulation (EU) 2017/745 (MDR), there is a documentation obligation for:

- Maintenance
- Incidents / near misses

12 Re-use

The points on **service and care & maintenance** specified in the operating instructions must be observed when cleaning the sliding board.

13 Waste disposal

13.1 Disposal of the device

The appliance and accessories, if present, should be disposed of in an environmentally friendly manner and in accordance with legal regulations. Please comply with the applicable **waste separation regulations!** If you are still unsure about this, please contact your local authority or waste disposal company.

13.2 Disposal of electrical components

*if electrical components are included in the medical device

In accordance with Directive 2012/19/EU - WEEE2, this medical device is classified as electrical equipment. All electrical components are free of non-authorized substances classified as harmful in accordance with RoHS II Directive 2011/65/EU. In addition, replaced electrical components must be disposed of in accordance with the European directives (see Directive 2012/19/EU - WEEE2).

13.3 Disposal of the packaging

The EU Waste Framework Directive 2008/98/EC is decisive for the handling of packaging disposal. Reusable materials must be recycled in accordance with national regulations.

14 EC Declaration of Conformity


As the manufacturer, we declare under our sole responsibility that our sliding boards comply with the essential requirements of the EC Directive for

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