Behandlungsliegen

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Instruction for use

Stand: 09/2022  
(Rev. 3.0)

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6.4	Operating the metal side rails .....	40
6.5	Operating the brake castors .....	41
6.6	Operating the central brake system .....	41
6.7	Operating the roller lifting system .....	41
6.8	Operating instructions .....	41
7	Ambient conditions .....	42
7.1	Storage conditions .....	42
7.2	Operating conditions .....	43
8	Technical data .....	44
9	Used materials .....	44
10	Service and care .....	45
11	Service life of the product .....	45
12	Disinfection .....	45
12.1	Specifications of detergents and disinfectants .....	45
13	Operational faults and solutions .....	46
14	Recommended accessories .....	46
15	Maintenance .....	47
15.1	Legal basis .....	47
15.2	Maintenance intervals .....	47
15.3	Spare parts .....	47
15.4	Notes on documentation .....	47
16	Disposal .....	48
16.1	Disposal of the device .....	48
16.2	Disposal of the electrical components .....	48
16.3	Disposal of the packaging .....	48
17	Declaration of Conformity .....	48

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## List of figures

Figure 1: Exemplary type plate ..... 32  
 Figure 2: Unscrewing the lying surface ..... 36  
 Figure 3: Lifting down the lying surface ..... 37  
 Figure 4: Screwing the lying surface ..... 37  
 Figure 5: 1-row hand control (1st row of keys: height adjustment) ..... 38  
 Figure 6: 2-row hand control (1st row of buttons: headrest adjustment, 2nd row of buttons: height adjustment) ..... 38  
 Figure 7: Operation of the step lever for hydraulic adjustment ..... 39  
 Figure 8: Fold-up side rails ..... 40  
 Figure 9: Metal side rails ..... 40  
 Figure 10: Exemplary castors in braked and unbraked condition ..... 41  
 Figure 11: Operation of the roller lifting system ..... 41

## List of tables

Table 1: Used symbols ..... 31  
 Table 2: Standards verification ..... 32  
 Table 3: Storage conditions ..... 42  
 Table 4: Operating conditions ..... 43  
 Table 5: Technical data ..... 44  
 Table 6: Operational faults and solutions ..... 46  
 Table 7: Recommended accessories ..... 46

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**Version: 09/2022  
(Rev. 3.0)**

## Table of contents

1	Foreword.....	6
2	General information.....	7
2.1	Symbols used .....	7
2.2	Type plate.....	9
2.3	Proof of standards.....	10
3	Safety instructions.....	11
4	General product description .....	12
4.1	Intended use .....	12
4.2	Indication .....	12
4.3	Contraindication.....	12
4.4	Equipment features.....	12
4.4.1	BL-001-0 BL-002-0 BL-003-0.....	12
4.4.2	BL-001-H BL-002-H BL-003-H.....	13
5	Assembly information .....	14
5.1	Basic information on installation .....	14
5.2	Disassembly/assembly instructions for the lying surface .....	14
5.3	Dismantling the treatment table.....	15
6	Service.....	16
6.1	Operation of the electrical movement.....	16
6.1.1	Automatic locking function .....	17
6.2	Operating the hydraulic movement.....	17
6.3	Operating the fold-up side rail .....	17
6.4	Operating the metal side rail .....	18
6.5	Operating the brake rollers.....	19
6.6	Operating the central braking system .....	19
6.7	Operating the roller lifting system .....	19
6.8	Operating instructions.....	20
7	Ambient conditions .....	21
7.1	Storage conditions.....	21
7.2	Operating conditions .....	21
8	Technical data .....	22
9	Materials used.....	23
10	Service and care.....	23
11	Useful life of the product .....	23

12	Disinfection .....	23
12.1	Specifications of cleaning agents and disinfectants: .....	23
13	Malfunctions and their rectification.....	24
14	Recommended accessories.....	24
15	Maintenance .....	25
15.1	Legal basis .....	25
15.2	Maintenance intervals.....	25
15.3	Spare parts .....	25
15.4	Notes on documentation .....	26
16	Waste disposal .....	26
16.1	Disposal of the device .....	26
16.2	Disposal of electrical components.....	26
16.3	Disposal of packaging.....	26
17	EC Declaration of Conformity .....	28

6.4	Operating the metal side rails .....	40
6.5	Operating the brake castors.....	41
6.6	Operating the central brake system .....	41
6.7	Operating the roller lifting system .....	41
6.8	Operating instructions .....	41
7	Ambient conditions .....	42
7.1	Storage conditions .....	42
7.2	Operating conditions.....	43
8	Technical data .....	44
9	Used materials .....	44
10	Service and care .....	45
11	Service life of the product.....	45
12	Disinfection .....	45
12.1	Specifications of detergents and disinfectants .....	45
13	Operational faults and solutions.....	46
14	Recommended accessories.....	46
15	Maintenance .....	47
15.1	Legal basis .....	47
15.2	Maintenance intervals .....	47
15.3	Spare parts.....	47
15.4	Notes on documentation .....	47
16	Disposal.....	48
16.1	Disposal of the device .....	48
16.2	Disposal of the electrical components .....	48
16.3	Disposal of the packaging.....	48
17	Declaration of Conformity.....	48

## List of illustrations

Figure 1: Exemplary type plate .....	9
Figure 2: Unscrewing the lying surface .....	14
Figure 3: Lifting down the lying surface .....	15
Figure 4: Fastening the lying surface .....	15
Figure 5: 1-row handset (1st row of buttons: height adjustment).....	16
Figure 6 : 2-row handset ( 1st row of buttons: head section adjustment, 2 n d row of buttons: height adjustment).....	16
Figure 7: Operating the pedal lever for hydraulic adjustment.....	17
Figure 8: Fold-up side rails .....	18
Figure 9: Metal side grille.....	18

Figure 10: Example castors in braked and unbraked condition.....19  
 Figure 11: Operating the roller lifting system.....19

## List of tables

Table 1: Symbols used .....9  
 Table 2: Verification of standards .....10  
 Table 3: Storage conditions.....21  
 Table 4: Operating conditions .....21  
 Table 5: Technical data .....22  
 Table 6: Malfunctions and their rectification.....24  
 Table 7: Recommended accessories .....24

## List of figures

Figure 1: Exemplary type plate .....32  
 Figure 2: Unscrewing the lying surface.....36  
 Figure 3: Lifting down the lying surface.....37  
 Figure 4: Screwing the lying surface .....37  
 Figure 5: 1-row hand control (1st row of keys: height adjustment) .....38  
 Figure 6: 2-row hand control (1st row of buttons: headrest adjustment, 2nd row of buttons: height adjustment) .....38  
 Figure 7: Operation of the step lever for hydraulic adjustment.....39  
 Figure 8: Fold-up side rails.....40  
 Figure 9: Metal side rails.....40  
 Figure 10: Exemplary castors in braked and unbraked condition.....41  
 Figure 11: Operation of the roller lifting system .....41

## List of tables

Table 1: Used symbols .....31  
 Table 2: Standards verification .....32  
 Table 3: Storage conditions .....42  
 Table 4: Operating conditions.....43  
 Table 5: Technical data .....44  
 Table 6: Operational faults and solutions.....46  
 Table 7: Recommended accessories .....46

## 1 Foreword

Dear customer!



We would like to thank you for the trust you have placed in us and for purchasing our product. We have manufactured this medical product with great care.

**Please read the operating instructions carefully before using the device for the first time and always keep them within easy reach.**

These instructions for use cannot cover all possible uses of the device. For further information or in the event of problems that are not described or not described in sufficient detail in these instructions for use, please contact your specialist dealer or medical supply retailer.

## 2 General information

### 2.1 Symbols used

	<p>This warning sign indicates all instructions that are important for safety. Non-compliance can lead to accidents or injuries.</p>
	<p>Manufacturer - Shows the manufacturer of the medical device in accordance with EU Directive 2017/745 The symbol, together with the name and address of the manufacturer (i.e. the person placing the medical device on the market), must be displayed in the immediate vicinity of the symbol. appear</p>
	<p>Conformity symbol according to 2017/745 of the Medical Devices Directive</p>
	<p>Medical device - Indicates the medical device provided by the manufacturer in accordance with EU Directive 2017/745</p>
	<p>Device type B according to IEC 601-1 (Special protection against electric shock)</p>
	<p>Protection class II device, protective insulation</p>
	<p>Dispose of electrical components in accordance with legal regulations. Do not dispose of in household waste!</p>
	<p>Date of manufacture - shows the date on which the medical device was manufactured.</p>
	<p>Article number - displays the manufacturer's article number so that the medical device can be identified.</p>

	<p>Serial number - displays the manufacturer's serial number so that a specific medical device can be identified.</p>
	<p>Sales partner - indicates the company that sells the medical device locally.</p>
	<p>Temperature limit - indicates the temperature limits to which the medical device can be safely exposed.</p>
	<p>Humidity, limit - indicates the humidity range to which the medical device can be safely exposed.</p>
	<p>Air pressure, limitation - indicates the range of air pressure to which the medical device can be safely exposed.</p>
	<p>Instructions for use Observe or electronic instructions for use - indicates to the user that it is necessary to follow the instructions for use.</p>
	<p>Unique identifier of a medical device - displays a carrier that contains information on a unique identifier of a medical device.</p>
	<p>Safe workload</p>
	<p>The product must be disposed of separately in the European Union. Disposal with normal household waste is not permitted.</p>

	<p>Minimum body measurements/weights of the patient</p>
--	---

Table 1: Symbols used

## 2.2 Type plate

The type plate is attached to the head end of the trolley frame. The product can be clearly identified by the type plate.


Figure 1: Exemplary type plate

Figure 1 shows an example nameplate. Please refer to the attached type plate for the exact specifications of your product

### 2.3 Proof of standards

The following national and international standards are used for the design and verification of the product, labeling and instructions for use.

Standard	Title	Issue
DIN EN 60601-1-6	Specification of suitability for use	2010
DIN EN 60601-1-2	Electromagnetic compatibility	2015
DIN EN ISO 10993	Biological evaluation of medical devices - Part 1: Assessment and verification	2010
DIN EN 1041	Provision of information by the manufacturer of a medical device	2008
DIN EN ISO 14971	Medical devices - Application of risk management to medical devices	2020

*Table 2: Verification of standards*

### 3 Safety instructions

- Before using the product, you should read these instructions for use carefully (see Medical Devices Operator Ordinance MPBetreibV §2). It contains important information for the safe and reliable use of the device. Keep the instructions for use in a safe place for future reference.
- Safety, reliability and performance are guaranteed if the following instructions are observed and the device is used properly. As the operator, you must comply with the Medical Devices Operator Ordinance (MPBetreibV).
- Ensure that children only have access to the treatment table under supervision and that no children are present in the area under the treatment table during operation.
- The table should only be set up by authorized personnel.
- The fuse protection on the installation side must not exceed 16A. Before connecting the appliance, please ensure that the voltage and frequency of your power supply correspond to the specifications on the rating plate.
- When selecting a location for the treatment table, make sure that the surface is level.
- Provide a suitable floor covering if the treatment table has to be moved frequently. Carpets, rugs and loose floor coverings can be damaged or make it difficult to move the treatment table.
- Connect the mains plug firmly to the mains socket. Lay the mains connection cable on the floor. Make sure that the treatment table (especially when moving it) does not rest on the cable with its castors. The cable must not be routed through the mechanism of the table base! (Danger of crushing)
- When the handsets are not in use, ensure that they are hanging on the treatment table and are not placed on the treatment table in order to avoid operating errors that could cause damage.
- If the patient is unattended, ensure that the treatment table is set at its lowest height to make it as easy as possible to get on and off.

## 4 General product description

### 4.1 Intended use

The treatment tables in the BL-001-0/BL-002-0/BL-003-0 series are equipped with electric height adjustment. They are operated using hand or foot switches. The drives for the adjustment functions consist of electromechanical linear motors with a maintenance-free

Permanent lubrication. The drives are operated via an air switch, which works without electric current but is actuated pneumatically. Treatment tables in the BL-001-H/BL-002-H/BL-003-H series are equipped with a hydraulic pump for adjusting the height of the table surface. This is operated via two foot levers on the chassis.

### 4.2 Indication

Treatment tables are aids whose use is indicated:

- for therapeutic treatments on patients who require a height adjustment of the lying surface to improve the therapy

### 4.3 Contraindication

The following patients are not permitted to use treatment tables:

- For severe cardiovascular problems
- In the event of severe dizziness, which no longer allows safe lying on the lying surface (e.g. falling down)
- For severe anxiety

### 4.4 Equipment features

#### 4.4.1 BL-001-0 BL-002-0 BL-003-0

The treatment tables have the following electrical functions:

- Electric height adjustment
- Possibility of locking the electrical functions

#### **Optional:**

- Electric headboard adjustment
- Electric leg adjustment

The drives for the adjustment functions consist of electromechanical linear motors with maintenance-free permanent lubrication. The drives are operated using a handset connected to the control unit via a spiral cable. The table has four single-brake castors.

The drives and the manual switch are electrically isolated from the mains voltage and are operated with a low voltage (DC 24 V).

#### **Optional:**

- Central braking system
- Roller lifting system
- Metal side rails (one or two-sided)
- Fold-up side rails (complete or head-end surround)
- Paper roll holder
- Nose slot in the upholstered panel

- Premium swivel castors
- Electric headboard adjustment
- Electric leg adjustment
- Headboard adjustment via a gas spring

#### 4.4.2 BL-001-H BL-002-H BL-003-H

The ...-H version does away with all electronic components. The height adjustment is realized via a hydraulic pump. It is operated via foot levers mounted on the left/right of the chassis.

**Optional:**

- Central braking system
- Roller lifting system
- Metal side rails (one or two-sided)
- Fold-up side rails (complete or head-end surround)
- Paper roll holder
- Nose slit in the upholstered panel
- Premium swivel castors
- Headboard adjustment via a gas spring

## 5 Assembly information

### 5.1 Basic information on installation

The treatment table should only be set up by authorized personnel. The fuse protection on the installation side must not exceed 16A.

Before connecting the appliance, please ensure that the voltage and frequency of your power supply correspond to the specifications on the rating plate.

When selecting a location for the treatment table, make sure that the surface is level. Provide a suitable floor covering if the treatment table has to be moved frequently. Carpets, rugs and loose floor coverings can be damaged or make it difficult to move.

Connect the mains plug firmly to the mains socket. Lay the mains connection cable on the floor. Make sure that the treatment table (especially when moving it) does not rest on the cable with its castors. The cable must not be routed through the mechanism of the table base! (Danger of crushing)

**Damage to the electrical mains cable caused by driving over or clamping can have fatal consequences.**

**Before moving the treatment table or dismantling it for transportation, the mains connection cable must be wound up and attached to the device provided on the chassis.**

### 5.2 Disassembly/assembly instructions for the lying surface

Estimated working time: 15 minutes

1. Unscrew the existing plate (six screws)

*Figure 2: Unscrewing the lying surface*

2. Place the new plate on the scissor mechanism

Figure 3: Lifting down the lying surface

3. Carefully screw the bolts diagonally into the thread of the new foam board.
  - a. Screw all bolts into the plate with your finger
  - b. Fasten the screws diagonally with the screwdriver

Figure 4: Fastening the lying surface

In the case of spare parts deliveries of lying surfaces, there may be deviations in the screw-on points due to the dimensional tolerance, which is why the lying surface is connected to the lying surface frame using conventional wood screws. These are included in the scope of delivery.

### 5.3 Dismantling the treatment table

If necessary, e.g. for transportation, the treatment tables can be dismantled with little effort as described but in reverse order. After dismantling the treatment table, it should be set up again by authorized personnel.

## 6 Operation

### 6.1 Operation of the electrical movement

(Electric adjustment for treatment tables in the BL-001-0/BL-002-0/BL-003-0 series)

Each upward or downward movement of the height adjustment is effected by pressing the respective button on the handset or foot switch. The buttons are marked with symbols according to their function.

Figure 5: 1-row handset (1st row of buttons: height adjustment)

**Do not exceed a maximum duty cycle of 6 minutes per hour!**

**If the table is not in use, it must be disconnected from the mains to prevent the height drives from moving in any case! This also applies to cleaning staff. They must be instructed accordingly when cleaning the table!**

If the treatment table has an additional electric head section adjustment, this must be operated using the corresponding handset in Figure 6.

Figure 6: 2-row handset (1st row of buttons: head section adjustment, 2nd row of buttons: height adjustment)

The electric height adjustment can also be operated using a foot switch. With (+), the bed will move upwards and with (-), it will move downwards. In this case, the automatic locking function must also be deactivated.

### 6.1.1 Automatic locking function


The treatment table has an automatic locking function ("Auto-Lock") for the electrical adjustment. When not in use, the electric adjustment of the table is always locked, which means that it cannot be operated using the hand control or foot control. The electric adjustment can be activated using a button combination on the switch. First press the "Up" button (1) and then the "Down" button (2). A short beep signals that the recliner has been unlocked and is now ready for use. After unlocking, the table is activated for 5 seconds. After this period, the table is automatically locked again. Locking is also signaled by a short beep. Unlocking is always carried out using the height button row.

## 6.2 Operation of the hydraulic movement

(Hydraulic adjustment for treatment tables in the BL-001-H/BL-002-H/BL-003-H series)

The lying surface can be adjusted in height by operating the pump levers on the left and right of the chassis.

When the right pump lever is pressed down, the lying surface is raised. Pressing down on the left foot pedal lowers the cushion plate.

*Figure 7: Operating the pedal lever for hydraulic adjustment*

## 6.3 Operation of the fold-up side rail

The folding side rail can be used to create an additional border around the lying surface. However, this border does not serve as protection against the patient falling out.

The folding process can be initiated by releasing the catch. To do this, pull out the black catch and turn it slightly so that the bolt no longer holds the side rail. Repeat this process for each side. The side rail can then be folded down. To return the side rail to an upright position, proceed in reverse order.

Figure 8: Fold-up side rails

**Persons must not sit on the side rails under any circumstances!**

#### 6.4 Operation of the metal side rail

Side rails can be used to stabilize independent patients on the relatively narrow beds or to support positioning cushions, for example. Side rails can also be used to prevent any risk of falling when repositioning or during treatment.

The position of the castor plays a decisive role when lowering the side rails. Please ensure that the metal side rails do not rest on the castors (this does not apply to the standard bed - dimensions 80x200 cm).

**Proper functioning of the latching fittings is a prerequisite; the proper latching of the fitting must also be monitored by the nursing staff!**

**Side rails on treatment tables are not suitable for restricting the mobility of persons requiring treatment.**

**Children in particular will try to overcome or bypass these barriers if there are only bars on the long sides of the loungers**

**In any case, patients must not be left on couches without competent supervision.**

**The treatment table may only be moved with the metal side rails raised.**

Figure 9: Metal side grille

## 6.5 Operating the brake rollers

A table in the standing version has four swivel castors with individual brakes. The table should always be braked at the installation site or during treatment using the castor lock. (see Figure 10)


Figure 10: Exemplary rollers in braked and unbraked condition

**Never leave the treatment table unrestrained in the room!**

## 6.6 Operating the central braking system

The treatment table should always be braked at the installation site by actuating the foot pedal of the central braking system so that the appliance is completely stable. In addition to the conventional implementation of the castors, there is also a solution using a central braking system. This can be purchased as an accessory. With the central brake system, all four castors are braked or released by activating a brake lever. The orientation of the brake lever indicates the status of the brake system. When the brake lever is in the horizontal position, the stretcher is unbraked. If the brake lever points in a direction rotated by 30°, the stretcher is braked.

## 6.7 Operation of the roller lifting system

When the foot lever is operated (foot lever pointing downwards, Figure 11 right), the castors are completely lifted off the floor, which means that the bed is completely stable and firm on its feet. The table cannot be moved in this position. To place the bed back on the castors, release the foot lever (foot lever is pointing upwards, Figure 11 left). The foot lever should not be operated in the lowest lying surface position, as it is more difficult to reach the lever in this position.


Figure 11: Operating the roller lifting system

## 6.8 Operating instructions

- After assembly and before the treatment table is used by a patient, check that all connections and the entire table are secure.
- Make sure that the drives are working properly.
- A treatment table that is not fully functional must be withdrawn from use immediately.
- Make sure that there are no objects, such as a wastepaper basket, side table, chair, etc., in the treatment table's movement area. side table, chair, etc. in the movement area of the treatment table.
- To avoid the risk of injury, no parts of the patient's body should protrude above the mattress base or feet should rest on the mattress base when the adjustment functions are being operated.
- Before moving the treatment table, the mains plug must be disconnected from the socket to prevent damage to the electrical system.
- The treatment table may only be moved with patients over thresholds with a maximum height of 2 cm.
- Ensure that the duty cycle is observed. Therefore, never make long and unnecessary electrical adjustments. If the thermal fuse in the control unit is triggered once after 6 min/h, the control unit must be replaced by an authorized specialist!

**The installation of additional devices such as insulin pumps, ventilators etc. is prohibited if no potential equalization has been created beforehand.**

**Any cables from additional devices must not be routed through the mechanism of the table base! (Danger of crushing)**

**If the table is not in use, it must be disconnected from the mains to prevent the height drives from moving in any case! This also applies to cleaning staff. They must be instructed accordingly when cleaning the table!**

## 7 Ambient conditions

### 7.1 Storage conditions



Storage temperature	min. -10 °C max. +50 °C	
Relative humidity	min. 20 % max. 80 %	
Air pressure (at altitude ≤ 3000 m)	min. 700 hPa max. 1060 hPa	

Table 3: Storage conditions

### 7.2 Operating conditions



Operating temperature	min. +5 °C max. +40 °C	
Relative humidity	min. 20 % max. 80 %	
Air pressure (at altitude ≤ 3000 m)	min. 700 hPa max. 1060 hPa	

Table 4: Operating conditions

## 8 Technical data

(Subject to change!)

Designation	BL-001(2)-H	BL-003-H	BL-001(2)-0	BL-003-0	BL-003-XS	BL-002-102	BL-003-S
Nominal voltage	230V~/ 50/60Hz						
Rated power	192 VA						
Device type B according to IEC 601-1							
Protection class							
IP protection class for drive components							
Control unit	IPX4						
Drive	IPX4						
Hand switch	IPX4						
Foot switch	IP X6						
Sound power level	50 dB(A)						
Duty cycle ED 10%	maximum 6 minutes/hour						
Height adjustment	49cm to 94 cm	49 cm-94	44 cm up to 95 cm	45 cm- 83 cm		44 cm to 95	
Safe workload	210 kg	260 kg	210 kg	300 kg		210 kg	260 kg
Dimensions of the lying surface (maximum)	90 cm x 200 cm	120 cm x 220 cm	90 cm x 200 cm	120 cm x 220 cm	110 x 200 cm	80 cm x 200 cm	120 cm x 220 cm
Total mass	appro x. 87 kg	appr ox. 96 kg	appro x. 65 kg	appro x. 95 kg	Approx. 90 kg	appro x. 75 kg	appro x. 95 kg
Height adjustment Stroke length	40 cm	40 cm	41 cm	38 cm	38 cm	41 cm	41 cm

Table 5: Technical data

**Repairs may only be carried out by ISKO specialist personnel or by persons authorized and trained by ISKO with comprehensive product knowledge. Any warranty and liability claims will be rejected if this provision is not complied with.**

## 9 Materials used

The treatment table is manufactured as a welded tubular steel construction. The surfaces are powder-coated or galvanized. The surfaces of this product are harmless to the skin from a health point of view.

## 10 Service and care

All household cleaners without ammonia and abrasive cleaners are permitted for cleaning the tube parts, the lying surface and the wooden parts with a damp cloth. Solvents (e.g. nitro) destroy the coating of the tubes and the lying surface!

Mechanical cleaning (e.g. scraping, sanding) or blast cleaning of the table is not permitted. All pivot points of the moving parts, including the bearing eyes on the adjuster, are fitted with maintenance-free plain bearings and must not be oiled or greased.

## 11 Useful life of the product

The expected service life is 7 years. Lack of maintenance and excessive strain on the product can significantly reduce the service life.

## 12 Disinfection

- To ensure that the examination table functions properly, every ISKO table should be cleaned, disinfected and checked after each use so that it can be used again immediately.
- Improper cleaning/disinfection of the table can be dangerous.
- Depending on the degree of soiling, we recommend cleaning the lounge with a damp cloth or similar.
- Disinfectants can be used for wipe and spray disinfection in their intended concentration. (see manufacturer's instructions)
- The dilution ratio recommended by the manufacturer in the respective instructions for use must be used.

**Solvents are not permitted.**

**Abrasives or scouring pads must not be used.**

### 12.1 Specifications of the cleaning agents and disinfectants:

- The working solutions should normally be used freshly prepared.
- Do not exceed or fall below the specified concentrations.
- They must not contain any corrosive or caustic components.
- They must not contain any substances that change the surface structure or the adhesion properties of the materials.
- Lubricants must not be attacked by the cleaning agent and disinfectant.

**Under no circumstances should soap or detergents be added to the disinfectant. There is a risk of explosion and fire if alcohol-based products are used over large areas.**

The use of unsuitable detergents and disinfectants may cause damage to the surface coating, for which ISKO KOCH GmbH cannot be held liable.

### 13 Malfunctions and their rectification

	Error	Measure
1	Motor switches off automatically and shows no reaction when the switch is pressed	Operating time too long, wait 20 minutes until the drive has cooled down
2	One motor or none of the motors reacts to the Manual switch actuation	Check the plug connection of the mains connection cable. Check the plug connection between the handset and control box. Check plug connection between motor cables and control box.
3	Beep from the control box → An actuator has reached its position lost	Reset/initialization Press and hold the row of height adjustment buttons simultaneously (really simultaneously) until the height drive moves completely. has initialized. The height drive retracts during initialization.

Table 6: Malfunctions and their rectification

Modifications, readjustments and repairs to the table that cannot be carried out in accordance with the above instructions may only be carried out by the manufacturer directly or by a workshop authorized by the manufacturer.

### 14 Recommended accessories

Article description	Order number
Foot switch lying on the floor (1 piece)	BL-101-0
Foot switch mounted on both sides of the chassis	BL-102-0
Roller lifting system	BL-501-0
1 metal side rail 145 cm, chrome-plated incl. brackets	BL-202-0
Metal side rail 145 cm, chrome-plated incl. brackets (1 pair)	BL-201-0
Paper roll holder for roll diameter max. 40 cm (supplied without paper roll)	BL-151-0
Nose slit in the upholstered panel incl. cover in the upholstered panel cover color	BL-306-1
Headboard adjustment with gas spring in place of the latching adjusters	BL-311-0
Electric head section adjustment with linear drive in place of the latching adjusters	BL-312-0
Nose slit in the head section incl. cover	BL-306-0
Mechanical foot section adjustment with detent adjusters (up to approx. +30°)	BL-999-3
Premium swivel castor (4 pieces)	BL-002-1W
Headboard adjustment with latching adjuster (without nose slot)	BL-310-1
Electric headboard adjustment with linear drive (without nose slit)	BL-312-3

Table 7: Recommended accessories

## 15 Maintenance

### 15.1 Legal basis

The Medical Device Regulation (EU) 2017/745 (MDR) and national laws and regulations require operators of medical devices to ensure that the medical device is in a safe operating condition throughout its entire service life.

### 15.2 Maintenance intervals

As a requirement of the Medical Devices Operator Ordinance §4 (Maintenance), a thorough visual inspection (1), a functional test (2) and a functional test (3) must be carried out after the medical device has been in operation for two years at the latest.

(2) and a current leakage test (3) in accordance with DIN EN 62353:2015-10.

(1) Special attention must be paid to the following points during the visual inspection:

- Tight fit of all screw connections
- Mobility of the pivot points
- Check the mains supply cable for crushing or shearing points
- Checking the strain relief of the mains supply cable

(2) Special attention must be paid to the following points during the functional test:

- Function of all electrically operated movements
- Fully extend and retract all the motors of the recliner until they switch off automatically. (Limit switches in the motors must switch off with an audible click)
- Functionality of the brakes
- Movability and function of the side rails
- Mobility of the triggers
- Checking the manual switch

**Functional tests and current leakage tests may only be carried out by ISKO specialist personnel or by persons authorized and trained by ISKO with comprehensive product knowledge.**

### 15.3 Spare parts

All spare parts for this medical device are available from ISKO KOCH GmbH, stating the serial number, order number and article number (these can be found on the type plate attached to the medical device).

To ensure that functional reliability and any warranty claims remain valid, only original ISKO KOCH GmbH parts should be used for spare parts.

ISKO KOCH GmbH

Egerländer Street 28

95448 Bayreuth

Phone: +49(0)921/150845-0 (Monday - Thursday 8:00 - 17:00 & Friday 8:00 - 16:00) Fax:

+49(0)921/150845-45

E-mail: [info@isko-koch.de](mailto:info@isko-koch.de)

## 15.4 Notes on documentation

According to the Medical Device Operator Ordinance and Medical Device Regulation (EU) 2017/745 (MDR), there is a documentation obligation for:

- Maintenance
- Incidents / near misses

## 16 Waste disposal

### 16.1 Disposal of the device

The appliance and accessories, if present, should be disposed of in an environmentally friendly manner and in accordance with legal regulations. Please comply with the applicable **waste separation regulations!** If you are still unsure about this, please contact your local authority or waste disposal company.

### 16.2 Disposal of electrical components

\*if electrical components are included in the medical device

In accordance with Directive 2012/19/EU - WEEE2, this medical device is classified as electrical equipment. All electrical components are free of non-authorized substances classified as harmful in accordance with RoHS II Directive 2011/65/EU. In addition, replaced electrical components must be disposed of in accordance with the European directives (see Directive 2012/19/EU - WEEE2).

### 16.3 Disposal of the packaging

The EU Waste Framework Directive 2008/98/EC is decisive for the handling of packaging disposal. Reusable materials must be recycled in accordance with national regulations.

## 17 EC Declaration of Conformity

As the manufacturer, we declare under our sole responsibility that our treatment tables comply with the essential requirements of the EC Directive for

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