Behandlungsliegen

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**Stand: 11/2022
(Rev. 2.1)**

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Status: 11/2022
(Rev. 2.1)

Table of contents

1	Foreword.....	4
2	General information.....	5
2.1	Symbols used	5
2.2	Nameplate.....	7
2.3	Proof of standards.....	8
3	Safety instructions.....	9
4	General product description	10
4.1	Intended use	10
4.2	Indication	10
4.3	Contraindication.....	10
4.4	Equipment features.....	10
5	Assembly information	11
5.1	Basic information on installation	11
5.2	Dismantling the standing board	11
6	Operation	12
6.1	Operating the straps	12
6.2	Operating the functions	13
6.3	Automatic locking function	13
6.4	Operating the emergency release.....	14
6.5	Operating the brake rollers.....	14
6.6	Operating instructions for special equipment and accessories	15
6.6.1	Electric headboard adjustment.....	15
6.6.2	Work and therapy table	15
6.6.3	Foot plates.....	16
6.6.4	Lying surface slot	17
6.6.5	Arm supports (support crutches)	17
6.7	Operating instructions.....	18
7	Ambient conditions	18
7.1	Storage conditions.....	18
7.2	Operating conditions	19
8	Technical data	19
9	Materials used.....	20
10	Service and care.....	20
11	Useful life of the product	20

12	Disinfection	20
12.1	Specifications of the cleaning agents and disinfectants:.....	20
13	Malfunctions and their rectification.....	21
14	Recommended accessories.....	21
15	Maintenance	22
15.1	Legal basis.....	22
15.2	Maintenance intervals.....	22
15.3	Spare parts.....	22
15.4	Notes on documentation	23
16	Waste disposal	23
16.1	Disposal of the device	23
16.2	Disposal of electrical components.....	23
16.3	Disposal of packaging.....	23
17	EC Declaration of Conformity	24

List of illustrations

Figure 1: Exemplary type plate	7
Figure 2: Procedure for looping in the webbing	12
Figure 3: 2-row handset (1st row of buttons: height adjustment, 2nd row of buttons: standing function)	13
Figure 4: Emergency release	14
Figure 5: Exemplary rollers in braked and unbraked condition.....	14
Figure 6: 3-row handset (1st row of buttons: height adjustment, 2nd row of buttons: standing function, 3rd row of buttons: head section adjustment).....	15
Figure 7: Worktop and tabletop	15
Figure 8: Mounting option for the work table.....	16
Figure 9: Position of the locking bolts of the foot plates.....	16
Figure 10: Adjusting cylinder for footplate adjustment.....	17
Figure 11: Lying surface slot.....	17

List of tables

Table 1: Symbols used	7
Table 2: Verification of standards	8
Table 3: Equipment features.....	10
Table 4: Storage conditions.....	18
Table 5: Operating conditions	19
Table 6: Technical data	19
Table 7: Malfunctions and their rectification.....	21
Table 8: Recommended accessories.....	22

1 Foreword

Dear customer!



We would like to thank you for the trust you have placed in us and for purchasing our product. We have manufactured this medical product with great care.


Please read the operating instructions carefully before using the device for the first time and always keep them within easy reach.

These instructions for use cannot cover all possible uses of the device. For further information or in the event of problems that are not described or not described in sufficient detail in these instructions for use, please contact your specialist dealer or medical supply retailer.

2 General information

2.1 Symbols used

	<p>This warning sign indicates all instructions that are important for safety. Non-compliance can lead to accidents or injuries.</p>
	<p>Manufacturer - Shows the manufacturer of the medical device in accordance with EU Directive 2017/745 The symbol, together with the name and address of the manufacturer (i.e. the person placing the medical device on the market), must be displayed in the immediate vicinity of the symbol. appear</p>
	<p>Conformity symbol according to 2017/745 of the Medical Devices Directive</p>
	<p>Medical device - Indicates the medical device provided by the manufacturer in accordance with EU Directive 2017/745</p>
	<p>Device type B according to IEC 601-1 (Special protection against electric shock)</p>
	<p>Protection class II device, protective insulation</p>
	<p>Dispose of electrical components in accordance with legal regulations. Do not dispose of in household waste!</p>
	<p>Date of manufacture - shows the date on which the medical device was manufactured.</p>
	<p>Article number - displays the manufacturer's article number so that the medical device can be identified.</p>

	<p>Serial number - displays the manufacturer's serial number so that a specific medical device can be identified.</p>
	<p>Sales partner - indicates the company that sells the medical device locally.</p>
	<p>Temperature limit - indicates the temperature limits to which the medical device can be safely exposed.</p>
	<p>Humidity, limit - indicates the humidity range to which the medical device can be safely exposed.</p>
	<p>Air pressure, limitation - indicates the range of air pressure to which the medical device can be safely exposed.</p>
	<p>Instructions for use Observe or electronic instructions for use - indicates to the user that it is necessary to follow the instructions for use.</p>
	<p>Unique identifier of a medical device - displays a carrier that contains information on a unique identifier of a medical device.</p>
	<p>Safe workload</p>
	<p>The product must be disposed of separately in the European Union. Disposal with normal household waste is not permitted.</p>

	<p>Minimum body measurements/weights of the patient</p>
--	---

Table 1: Symbols used

2.2 Type plate

The type plate is attached to the head end of the trolley frame. The type plate allows the product to be clearly identified.


Figure 1: Exemplary type plate

Figure 1 shows an example nameplate. Please refer to the attached type plate for the exact specifications of your product

2.3 Proof of standards

The following national and international standards are used for the design and verification of the product, labeling and instructions for use.

Standard	Title	Issue
DIN EN 60601-1-6	Specification of suitability for use	2010
DIN EN 60601-1-2	Electromagnetic compatibility	2015
DIN EN ISO 10993	Biological evaluation of medical devices - Part 1: Assessment and verification	2010
DIN EN 1041	Provision of information by the manufacturer of a medical device	2008
DIN EN ISO 14971	Medical devices - Application of risk management to medical devices	2020

Table 2: Verification of standards

3 Safety instructions

- Before using the product, you should read these instructions for use carefully (see Medical Devices Operator Ordinance MPBetreibV §2). It contains important information for the safe and reliable use of the device. Keep the instructions for use in a safe place for future reference.
- Safety, reliability and performance are guaranteed if the following instructions are observed and the device is used properly. As the operator, you must comply with the Medical Devices Operator Ordinance (MPBetreibV).
- The standing boards are suitable for both home care and clinical use, with a maximum patient weight of 135 kg.
- Ensure that children only have access to the standing board under supervision and that no children are in the danger zone under the standing board during operation.
- The standing board should only be set up by authorized personnel.
- The fuse protection on the installation side must not exceed 16A. Before connecting the appliance, please ensure that the voltage and frequency of your power supply correspond to the specifications on the rating plate.
- Make sure that the standing board is placed on a level surface.
- Provide a suitable floor covering if the standing board has to be moved frequently. Carpets, rugs and loosely laid floor coverings can be damaged or make it difficult to move the standing board.
- Connect the mains plug firmly to the mains socket. Lay the mains connection cable on the floor. Make sure that the standing board (especially when moving it) does not rest on the cable with its castors. The cable must not be routed through the mechanics of the standing board base! (Danger of crushing)
- When the handsets are not in use, make sure that they are hanging on the standing board and are not placed in the standing board to avoid operating errors that could cause damage.

The patient must be strapped in with all three belts before standing up. Only operate the stand-up function if an assistant is present. The assistant or the attending person must perform the operations and stand in front of the patient. It is essential to prevent any part of the body of any person getting into the rear part of the device (risk of crushing)!

4 General product description

4.1 Intended use

The ST202 standing board is a medical aid for patients with severe limitations in their ability to stand or walk. The ST202 standing board differs from its predecessor model ST101 in that it has an additional electric height adjustment.

The

Standing board was developed for the therapy of patients who need to get used to standing again, or for better urination and circulation training for bedridden patients.

4.2 Indication

People between 140 and 200 cm tall with various physical and/or psychological limitations can stand safely, firmly and without fatigue every day. Daily standing stabilizes the circulation, prevents bone decalcification (osteoporosis), corrects contractures of the hip and knee joints, stimulates bowel and bladder functions as well as breathing and kidney activity. Targeted promotion of head control is optimized by the precise adjustment options at all levels.

Standing boards are aids whose use is indicated:

- Paraplegia
- Paresis
- Spasticity
- Multiple disabilities of various kinds and hypertension
- For people who are severely or completely restricted in their ability to walk or stand

4.3 Contraindication

The following patients are not permitted to use standing boards:

- For pressure sores, especially in the hip and leg area
- In case of extreme deformation and non-weight-bearing capacity of the lower extremities
- For severe cardiovascular problems
- In the case of severe dizziness which makes it impossible to stand up independently or partially independently
- For severe anxiety

4.4 Equipment features

These are the electrical functions of the standing board:

Type	Electrical Height adjustment	Electrical Standing function	Electrical Headboard adjustment
SB-111-0	Yes	Yes	optional

Table 3: Equipment features

The drives for the adjustment functions consist of electromechanical linear motors with maintenance-free permanent lubrication. The drives are operated using a hand switch connected to the control unit via a spiral cable. The standing board has four single-brake castors.

The drives and the manual switch are electrically isolated from the mains voltage and are operated with a low voltage (DC 24 V).

5 Assembly information

5.1 Basic information on installation

The standing board should only be set up by authorized personnel. The fuse protection on the installation side must not exceed 16A.

Before connecting the appliance, please ensure that the voltage and frequency of your power supply correspond to the specifications on the rating plate.

When selecting a location for the standing board, make sure that the surface is level. Provide a suitable floor covering if the standing board has to be moved frequently. Carpets, rugs and loose floor coverings can be damaged or make it difficult to move.

Connect the mains plug firmly to the mains socket. Lay the mains connection cable on the floor. Make sure that the standing board (especially when moving it) does not rest on the cable with its castors. The cable must not be routed through the mechanics of the standing board base! (Danger of crushing)

Damage to the electrical mains cable caused by driving over or clamping can have fatal consequences.

Before moving the standing board or dismantling it for transportation, the mains connection cable must be wound up and attached to the device provided on the chassis.

5.2 Dismantling the standing board

If necessary, e.g. for transportation, the standing boards can be dismantled with little effort as described but in reverse order. Reassembly after dismantling the standing board should be carried out by authorized personnel.

6 Operation

6.1 Operating the straps

Before setting up the lying surface in the standing position, the patient must be secured around the chest, hips and knees with a retaining strap

When looping the belt through the belt tongue, please ensure that you loop the belt into the belt tongue from the correct side.

Figure 2: Procedure for looping in the webbing

Push the webbing under the belt tongue through slot 1. Pressing the belt tongue together creates a second opening (slot 2) through which the belt is guided backwards again. (see figure 2)

The length of the belt can thus be individually adjusted to each patient.

The patient must not be placed in the standing position until they have been adequately secured with the belt system!

6.2 Operating the functions

Each upward or downward movement of the standing function and the height adjustment of the lying surface is effected by pressing the respective buttons on the hand control marked with symbols. It should be noted that the standing function can only be initiated when the height adjustment has been moved all the way up. This is achieved using software that first moves the height adjustment to its top position when the "Up" standing function is activated.

If the standing board is in the standing position, the lying surface can be returned to the horizontal position by pressing 2 buttons. Here you can either press "Lower lying height" or "Standing function down". The end position when pressing "Standing function down" is a horizontal lying position in the highest position. Whereas the end position when pressing "Lying height lower" is a horizontal lying position in the lowest position.


Figure 3: 2-row handset (1st row of buttons: height adjustment, 2nd row of buttons: standing function)

The second row of buttons enables standing and tilting back from any position of the standing board. To start & tilt back this function, the "Stand up" & "Stand down" buttons must always be pressed twice in quick succession.

To start the stand-up function, press this button twice in quick succession.

When tilting the lying surface back, pay attention to the clamping points underneath the lying surface.

6.3 Automatic locking function

The standing board has an automatic locking function ("Auto-Lock") for the electric height adjustment. When not in use, the electric height adjustment of the standing board is always locked, which means that it cannot be operated using the handset. The electric adjustment can be activated using a button combination on the switch. First press the "Up" button (1) and then the "Down" button (2) of the corresponding height button row. A short beep indicates that the recliner has been unlocked and is now ready for use. After unlocking, the table is activated for 5 seconds. After this period, the table is automatically locked again.

Locking is also signaled by a short beep. Unlocking is always carried out via the height button row.

6.4 Operating the emergency release

In the event of a defect or a power failure, a manual emergency lowering system was installed in the standing board to bring the padded plate with the patient into a horizontal position.

A clamping lever (see Figure 4) is located in the center of the lifting motor mount. To move the lying surface by hand from the standing position to the horizontal position, proceed as follows:

- Pull out the safety spring plug on the motor bolt
- Hold the motor at the end of the push rod with one hand, then operate the clamping lever in the direction of the arrow with the other hand.
- **Please note:** The motor responsible for the standing function must be pressed down slightly so that it can be swung downwards.
- The lying surface can now be swiveled slowly (with the patient) into the horizontal position

Thanks to two built-in gas springs, the lying surface can be easily tilted back into the horizontal position with the patient without much effort. For safety reasons, at least two assistants should be present during manual lowering.


Figure 4: Emergency release

6.5 Operating the brake rollers

The bed has four swivel castors with individual brakes. The bed should always be braked at the installation site using the castor lock. (see Figure 5)


Figure 5: Exemplary rollers in braked and unbraked condition

6.6 Operating instructions for special equipment and accessories

6.6.1 Electric headboard adjustment

Art. No. SB-170-0

This accessory can be used to adjust the headboard electrically. The headboard has an angle adjustment range of -30° - $+60^{\circ}$.

Figure 6: 3-row handset (1st row of buttons: height adjustment, 2nd row of buttons: standing function, 3rd row of buttons: head section adjustment)

The head section adjustment is operated using a 3-row handset. The third row of buttons is used for headrest adjustment.

Before the standing board is tilted back into the horizontal position after standing, the electrically adjustable head section must be in the starting position, i.e. in line with the lying surface.

If the lying surface has reached an angle of approx. $+25^{\circ}$ when starting the standing process, the head section can no longer be positively adjusted. This is to prevent the patient from being moved to unfavorable positions by operating the head section.

6.6.2 Work and therapy table

Art. No. SB-141-0, SB-141-1

Two brackets are clamped to the frame at the desired height using the star grips.

Figure 7: Worktop and tabletop

The work table is then inserted with the chrome tubes into the holders, pushed towards the patient and secured using the star handles. (see Figure 8)

The table can be assembled by one person without much effort, as the brackets are fitted individually and in advance. We recommend placing the lying surface in the standing position.


Figure 8: Mounting option for the work table

6.6.3 Foot plates

Art. No. SB-121-0 (angle and length adjustable)

Length adjustment:

Release the 2 locking bolts (pos. 1, Figure 9) underneath the footplate by pulling them and bring them into the freewheel position by turning them slightly to the left or right. Make the desired adjustment by pulling the plate away from the lying surface or pushing it towards the lying surface. The bolts must now be returned to the locking slot by turning them. A slight tug on the footplate makes it easier to engage.

Figure 9: Position of the locking bolts of the foot plates

Please always ensure that the locking bolts are engaged.

Angle adjustment:

There is an adjusting cylinder between the footplate and the end of the lying surface (viewed from above). (Pos. 2, Figure 10) The footplate can be tilted to the desired angle by turning the adjusting cylinder to the left or right.

Figure 10: Adjusting cylinder for footplate adjustment

6.6.4 Lying surface slot

Art. No. SB-171-0 (incl. strap for individual fixation of the legs)

The special belt (see Figure 11) can be used to fix the patient's legs and feet in position before "standing". The length of the belt can be adjusted individually for each patient.

Figure 11: Lying surface slot

6.6.5 Armrests (support crutches)

Art. No. SB-121-0

Arm supports (support crutches) are bars that are attached to the lying surface and allow the patient to support themselves with both hands. The two holders of the arm supports are clamped to the star handles at the required height. The support crutches are then inserted into the holders and the star handles are tightened.

6.7 Operating instructions

- After assembly and before the table is used by a patient, check that all connections and the entire table are secure.
- To avoid the risk of injury, no parts of the patient's body may protrude above the mattress base or feet may rest on the mattress base while the adjustment functions are being operated.
- Make sure that the drives are working properly.
- A standing board that is not fully functional must be withdrawn from use immediately.
- Make sure that there are no objects, such as a waste paper basket, side table, chair, etc., in the movement area of the standing board. side table, chair, etc. in the movement area of the standing board.
- Before moving the standing board, the mains plug must be disconnected from the socket to prevent damage to the electrical system.
- Ensure that the duty cycle is observed. Therefore, never make long and unnecessary electrical adjustments. If the thermal fuse in the control unit is triggered once after 6 min/h, the control unit must be replaced by an authorized specialist!

The standing board is not intended for transporting patients!

The installation of additional devices such as insulin pumps, ventilators etc. is prohibited if no potential equalization has been created beforehand.

Possible cables from additional devices must not through the mechanics of the standing board substructure! (Danger of crushing)

7 Ambient conditions

7.1 Storage conditions

Storage temperature	min. -10 °C max. +50 °C	
Relative humidity	min. 20 % max. 80 %	
Air pressure (at altitude ≤ 3000 m)	min. 700 hPa max. 1060 hPa	

Table 4: Storage conditions

7.2 Operating conditions



Operating temperature	min. +5 °C max. +40 °C	
Relative humidity	min. 20 % max. 80 %	
Air pressure (at altitude ≤ 3000 m)	min. 700 hPa max. 1060 hPa	

Table 5: Operating conditions

8 Technical data

(Subject to change!)

	SB-111-0
Nominal voltage	~230V / 50Hz
Rated power	1035 VA
Duty cycle ED 10%	maximum 6 minutes/hour
Device type B according to IEC 601-1	
Protection class	
IP protection class	IPX4
Maximum patient weight	135 kg
Safe workload	170 kg
Total length	205cm
Total width	73 cm, 100 cm with attached therapy table
Dimensions of the upholstered panel	73 cm x 185 cm
Height adjustment	44 cm to 84 cm
Swivel range (standing function)	max. 85° (stepless adjustment)
Height during complete approached standing function	220cm

Table 6: Technical data

Repairs may only be carried out by ISKO specialist personnel or by persons authorized and trained by ISKO with comprehensive product knowledge. Any warranty and liability claims will be rejected if this provision is not complied with.

9 Materials used

The treatment table is manufactured as a welded tubular steel construction. The surfaces are powder-coated or galvanized. The surfaces of this product are harmless to the skin from a health point of view.

10 Service and care

All household cleaners without ammonia and abrasive cleaners are permitted for cleaning the tube parts, the lying surface and the wooden parts with a damp cloth. Solvents (e.g. nitro) destroy the coating of the tubes and the lying surface!

Mechanical cleaning (e.g. scraping, sanding) or blast cleaning of the table is not permitted. All pivot points of the moving parts, including the bearing eyes on the adjuster, are fitted with maintenance-free plain bearings and must not be oiled or greased.

11 Useful life of the product

The expected service life is 7 years. Lack of maintenance and excessive strain on the product can significantly reduce the service life.

12 Disinfection

- To ensure that the examination table functions properly, every ISKO table should be cleaned, disinfected and checked after each use so that it can be used again immediately.
- Improper cleaning/disinfection of the table can be dangerous.
- Depending on the degree of soiling, we recommend cleaning the lounge with a damp cloth or similar.
- Disinfectants can be used for wipe and spray disinfection in their intended concentration. (see manufacturer's instructions)
- The dilution ratio recommended by the manufacturer in the respective instructions for use must be used.

Solvents are not permitted.

Abrasives or scouring pads must not be used.

12.1 Specifications of the cleaning agents and disinfectants:

- The working solutions should normally be used freshly prepared.
- Do not exceed or fall below the specified concentrations.
- They must not contain any corrosive or caustic components.
- They must not contain any substances that change the surface structure or the adhesion properties of the materials.
- Lubricants must not be attacked by the cleaning agent and disinfectant.

Under no circumstances should soap or detergents be added to the disinfectant. There is a risk of explosion and fire if alcohol-based products are used over large areas.

The use of unsuitable detergents and disinfectants may cause damage to the surface coating, for which ISKO KOCH GmbH cannot be held liable.

13 Malfunctions and their rectification

	Error	Measure
1	Motor switches off automatically and shows no reaction when Actuating the switch	Operating time too long, wait 20 minutes until the drive has cooled down
2	One motor or none of the motors reacts to the Manual switch actuation	Check the plug connection of the mains connection cable. Check the plug connection between the handset and control box. Check plug connection between motor cables and control box.
3	Beep from the control box → An actuator has reached its position lost	Reset/initialization Press and hold the second row of buttons simultaneously (really simultaneously) until the interrupted signal tone changes to a continuous tone (after approx. 5 seconds). Immediately after the reset (simultaneous pressing of the 2nd row of buttons), the first row of buttons is held down simultaneously for initialization until a long signal tone sounds. In this process the motors may search for their end position.

Table 7: Malfunctions and their rectification

Modifications, readjustments and repairs to the table that cannot be carried out in accordance with the above instructions may only be carried out by the manufacturer directly or by a workshop authorized by the manufacturer.

14 Recommended accessories

Article description	Order number
Foot plates, individually angle and length adjustable (1 pair)	SB-121-0
Support crutches (1 pair)	SB-131-0
Wooden worktable, with body cut-out	SB-141-0
Plexiglas worktable, with body cut-out	SB-141-1
Electric headboard adjustment	SB-170-0
Lying surface slot incl. belt for securing the legs and feet	SB-171-0
Upholstery plate with nose slit, approx. 8 cm x 18 cm	SB-187-0
Upholstered panel with nose slit, approx. 8 cm x 30 cm	SB-187-2
Cover for nose slit, approx. 8 cm x 18 cm	SB-187-1

Cover for nose slit, approx. 8 cm x 30 cm	SB-187-3
Head fixation	SB-181-0
Raised chassis, to increase underride clearance	SB-188-0

Table 8: Recommended accessories

15 Maintenance

15.1 Legal basis

The Medical Device Regulation (EU) 2017/745 (MDR) and national laws and regulations require operators of medical devices to ensure that the medical device is in a safe operating condition throughout its entire service life.

15.2 Maintenance intervals

As a requirement of the Medical Devices Operator Ordinance §4 (Maintenance), a thorough visual inspection (1), a functional test (2) and a functional test (3) must be carried out after the medical device has been in operation for two years at the latest.

(2) and a current leakage test (3) in accordance with DIN EN 62353:2015-10.

- (1) Special attention must be paid to the following points during the visual inspection:
- Tight fit of all screw connections
 - Mobility of the pivot points
 - Check the mains supply cable for crushing or shearing points
 - Checking the strain relief of the mains supply cable
- (2) Special attention must be paid to the following points during the functional test:
- Function of all electrically operated movements
 - Fully extend and retract all motors on the standing board until they switch off automatically. (Limit switches in the motors must switch off with an audible click)
 - Functionality of the brakes
 - Movability and function of the side rails
 - Mobility of the triggers
 - Checking the manual switch

Functional tests and current leakage tests may only be carried out by ISKO specialist personnel or by persons authorized and trained by ISKO with comprehensive product knowledge.

15.3 Spare parts

All spare parts for this medical device are available from ISKO KOCH GmbH, quoting the serial number, order number and article number (these can be found on the type plate attached to the medical device).

To ensure that functional reliability and any warranty claims remain valid, only original ISKO KOCH GmbH parts should be used for spare parts.

ISKO KOCH GmbH

Egerländer Street 28

95448 Bayreuth

Phone: +49(0)921/150845-0 (Monday - Thursday 8:00 - 17:00 & Friday 8:00 - 16:00)

15.4 Notes on documentation

According to the Medical Device Operator Ordinance and Medical Device Regulation (EU) 2017/745 (MDR), there is a documentation obligation for:

- Maintenance
- Incidents / near misses

16 Waste disposal

16.1 Disposal of the device

The appliance and accessories, if present, should be disposed of in an environmentally friendly manner and in accordance with legal regulations. Please comply with the applicable **waste separation regulations!** If you are still unsure about this, please contact your local authority or waste disposal company.

16.2 Disposal of electrical components

*if electrical components are included in the medical device

In accordance with Directive 2012/19/EU - WEEE2, this medical device is classified as electrical equipment. All electrical components are free of non-authorized substances classified as harmful in accordance with RoHS II Directive 2011/65/EU. In addition, replaced electrical components must be disposed of in accordance with the European directives (see Directive 2012/19/EU - WEEE2).

16.3 Disposal of the packaging

The EU Waste Framework Directive 2008/98/EC is decisive for the handling of packaging disposal. Reusable materials must be recycled in accordance with national regulations.

17 EC Declaration of Conformity

As the manufacturer, we declare under our sole responsibility that our standing boards comply with the basic requirements of the EC Directive for

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