Levatrain (Max) Levaflex

ISK^omed[®]

Stehgeräte

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Stand: 09/2022
(Rev. 2.0)

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Status: 09/2022
(Rev. 2.0)

Table of contents

1	Foreword.....	5
2	General information.....	5
2.1	Symbols used	5
2.2	Nameplate.....	7
2.3	Proof of standards.....	8
3	Safety instructions.....	9
4	General product description	10
4.1	Intended use	10
4.2	Indication	10
4.3	Contraindication.....	10
4.4	Equipment features.....	10
4.5	Scope of delivery	10
5	Assembly information	11
5.1	Basic information on installation	11
6	Operation	11
6.1	Body size	11
6.2	Back pad or hip pad.....	12
6.3	Foot and heel support.....	13
6.4	Knee pads.....	13
6.5	Upper body pad.....	13
6.5.1	Adjusting the backrest cushion	13
6.5.2	Lateral armrest.....	14
6.5.3	Adjusting the width of the lateral hip supports	14
6.5.4	Attaching and adjusting the headrest.....	14
6.6	Breast pad	14
6.7	General instructions for use when standing up with the Levatrain	15
6.8	General instructions for use when standing up with the Levaflex.....	16
6.9	Behavior in the event of a power failure.....	18
6.10	Operating the brake rollers	18
7	Ambient conditions	19
7.1	Storage conditions.....	19
7.2	Operating conditions	19
8	Technical data subject to change!.....	19
9	Materials used.....	20

10	Service and care	20
11	Useful life of the product	20
12	Disinfection	20
12.1	Specifications of cleaning agents and disinfectants:	21
13	Malfunctions and their rectification.....	21
14	Recommended accessories	21
15	Maintenance	22
15.1	Legal basis	22
15.2	Maintenance intervals.....	22
15.3	Spare parts	22
15.4	Notes on documentation	23
16	Re-use	23
17	Waste disposal	23
17.1	Disposal of the device	23
17.2	Disposal of electrical components.....	23
17.3	Disposal of the packaging.....	24
18	EC Declaration of Conformity	24

List of illustrations

Figure 1: Exemplary type plate	7
Figure 2: Adjusting the table height.....	12
Figure 3: Clamping lever for table height adjustment.....	12
Figure 4: Position of the hip pad in the buttocks	13
Figure 5: Adjusting the heel switches and toe fixation	13
Figure 6: Adjusting the knee pads.....	13
Figure 7: Upper body pad.....	14
Figure 8: Breast pad	15
Figure 9: Patient restraint belt	15
Figure 10: Position of emergency stop pushbutton.....	16
Figure 11: Manual switch for Levatrain and Levatrain Max.....	16
Figure 12: Adjusting the castors.....	16
Figure 13: Manual switch for Levaflex	17
Figure 14: Dismantling process	17
Figure 15: Exemplary rollers in braked and unbraked condition	18

List of tables

Table 1: Symbols used	7
Table 2: Verification of standards	8
Table 3: Procedure for standing up	17
Table 4: Reset procedure.....	18
Table 5: Storage conditions.....	19
Table 6: Operating conditions	19
Table 7: Technical data	20
Table 8: Malfunctions and their rectification.....	21
Table 9: Recommended accessories.....	22

1 Foreword

Dear customer!



We would like to thank you for the trust you have placed in us and for purchasing our product. We have manufactured this medical device with great care.


Please read the operating instructions carefully before using the device for the first time and always keep them within easy reach.

These instructions for use cannot cover all possible uses of the device. For further information or in the event of problems that are not described or not described in sufficient detail in these instructions for use, please contact your specialist dealer or medical supply retailer.

2 General information

2.1 Symbols used

	<p>This warning sign indicates all instructions that are important for safety. Non-compliance can lead to accidents or injuries.</p>
	<p>Manufacturer - Shows the manufacturer of the medical device in accordance with EU Directive 2017/745 The symbol, together with the name and address of the manufacturer (i.e. the person placing the medical device on the market), must be displayed in the immediate vicinity of the symbol. appear</p>
	<p>Conformity symbol according to 2017/745 of the Medical Devices Directive</p>
	<p>Medical device - Indicates the medical device provided by the manufacturer in accordance with EU Directive 2017/745</p>
	<p>Device type B according to IEC 601-1 (Special protection against electric shock)</p>
	<p>Protection class II device, protective insulation</p>

	<p>Dispose of electrical components in accordance with legal regulations. Do not dispose of in household waste!</p>
	<p>Date of manufacture - shows the date on which the medical device was manufactured.</p>
	<p>Article number - displays the manufacturer's article number so that the medical device can be identified.</p>
	<p>Serial number - displays the manufacturer's serial number so that a specific medical device can be identified.</p>
	<p>Sales partner - indicates the company that sells the medical device locally.</p>
	<p>Temperature limit - indicates the temperature limits to which the medical device can be safely exposed.</p>
	<p>Humidity, limitation - displays the humidity range range, to which the medical device can be safely exposed to.</p>
	<p>Air pressure, limitation - indicates the range of air pressure to which the medical device can be safely exposed.</p>
	<p>Instructions for use Observe or electronic instructions for use - indicates to the user that it is necessary to follow the instructions for use.</p>

	<p>Unique identifier of a medical device - displays a carrier that contains information on a unique identifier of a medical device.</p>
	<p>Safe workload</p>
	<p>Permissible patient weight</p>
	<p>The product must be disposed of separately in the European Union. Disposal with normal household waste is not permitted.</p>
	<p>Minimum body measurements/weights of the patient</p>

Table 1: Symbols used

2.2 Type plate

The type plate is attached to the head end of the trolley frame. The product can be clearly identified by the type plate.

Figure 1: Exemplary type plate

Figure 1 shows an example nameplate. Please refer to the attached type plate for the exact specifications of your product

2.3 Proof of standards

The following national and international standards are used for the design and verification of the product, labeling and instructions for use.

Standard	Title	Issue
DIN EN 60601-1-6	Specification of suitability for use	2010
DIN EN 60601-1-2	Electromagnetic compatibility	2015
DIN EN ISO 10993	Biological evaluation of medical devices - Part 1: Assessment and verification	2010
DIN EN 1041	Provision of information by the manufacturer of a medical device	2008
DIN EN ISO 14971	Medical devices - Application of risk management to medical devices	2020

Table 2: Verification of standards

3 Safety instructions

- Before using the standing device, you should read these instructions for use carefully (see Medical Devices Operator Ordinance MPBetreibV §2). It contains important information for the safe and reliable use of the device. Keep the instructions for use in a safe place for future reference.
- Safety, reliability and performance are guaranteed if the following instructions are observed and the device is used properly. As the operator, you must comply with the Medical Devices Operator Ordinance (MPBetreibV).
- The device may only be operated by medically trained and instructed personnel at the start of upright training. Operators must be trained in the use of the device and their training and level of knowledge must enable them to assess the risks for the patient. The trained personnel are responsible for the assessment for subsequent independent use by the user (patient), whereby an assistant must remain within call or sight.
- All clamp and screw connections must be checked regularly
- During use, 1 or 2 assistants must be available to support the user.
- Only use aids on a firm, level and horizontal surface
- Do not use the aid outdoors. The device may only be used in protected domestic and clinical areas.
- Leaning out far poses a risk of falling
- **The aid must not be used as a means of transportation**
- Only park the aid with the wheels locked
- In the event of a fault, the standing appliance must be withdrawn from use immediately
- Parts of the standing appliance can catch fire if handled improperly
- When standing, the user must always wear sturdy shoes; barefoot use is not permitted. Unless expressly instructed by the therapist
- Ensure that children only have access to the standing appliance under supervision.
- The standing device should only be set up by authorized personnel.
- The fuse protection on the installation side must not exceed 16A. Before connecting the appliance, please ensure that the voltage and frequency of your power supply correspond to the specifications on the rating plate.
- When choosing a location, make sure the surface is level
- Provide a suitable floor covering if the standing appliance has to be moved frequently. Carpets, rugs and loose floor coverings can be damaged or make it difficult to move the appliance.
- Connect the mains plug firmly to the mains socket. Lay the mains connection cable on the floor. Make sure that the upright appliance (especially when moving it) does not stand on the cable with its castors. The cable must not be routed through the mechanics of the substructure! (Danger of crushing)
- Damaged mains cables can lead to life-threatening situations. They must be replaced immediately.
- Check the mains cable for damage at regular intervals (weekly).
- Ensure that the electrical specifications of the appliance match the local conditions at the installation site.
- When the handsets are not in use, ensure that they are attached to the standing appliance to prevent incorrect operation that could cause damage.

4 General product description

4.1 Intended use

The Levatrain standing device is used exclusively for standing training for people who are unable to walk or have difficulty walking, for independent use with an assistant. Standing training is used in inpatient and outpatient care.

rehabilitation, as well as in the home environment. The design of our standing devices enables independent, stable and safe standing. An electric belt retractor brings the patient out of the wheelchair and into a secure standing position.

4.2 Indication

People between 140 and 200 cm tall with various physical and/or mental limitations can stand safely, firmly and without fatigue every day. Daily standing stabilizes the circulation, prevents bone decalcification (osteoporosis), corrects contractures of the hip and knee joints, stimulates bowel and bladder functions as well as breathing and kidney activity. The targeted promotion of head control is optimized by the precise adjustment options at all levels.

4.3 Contraindication

The following patients are not permitted to use the standing device:

- For pressure sores, especially in the hip and leg area
- In case of extreme deformation and non-weight-bearing capacity of the lower extremities
- For severe cardiovascular problems
- In the case of severe dizziness which makes it impossible to stand up independently or partially independently
- For severe anxiety

4.4 Equipment features

The Levatrain standing device has the following features:

- Table cut-out with ergonomic belly pad
- Foot plate with heel cups and toe fixation
- Double swivel castors with total locking
- Table height adjustment with gas spring support
- Knee pads adjustable in height and depth, as well as with horizontal and angle adjustment
- Stainless steel handrails height-adjustable in parallel
- Belt retractor system available electrically
- Simple maintenance and repair of the electrical or mechanical belt retractor system

4.5 Scope of delivery

- Pre-installed stable base frame
- Support table with height adjustment and electric automatic belt retractor
- Height and depth-adjustable back cushions and knee pads
- adjustable heel cups with toe Velcro straps

5 Assembly information

5.1 Basic information on installation

The standing appliance should only be set up by authorized personnel. The fuse protection on the installation side must not exceed 16A.

Before connecting the appliance, please ensure that the voltage and frequency of your power supply correspond to the specifications on the rating plate.

Ensure that the floor is level when selecting a location for the standing device. Provide a suitable floor covering if the standing device has to be moved frequently. Carpets, rugs and loose floor coverings can be damaged or make it difficult to move.

Connect the mains plug firmly to the mains socket. Lay the mains connection cable on the floor. Make sure that the standing appliance (especially when moving it) does not rest on the cable with its castors.

Damage to the electrical mains cable caused by driving over or clamping can have fatal consequences.

6 Operation

All presettings may only be made on the empty standing device without the patient standing in it. The settings must then be checked again with the patient standing in the aid and corrected if necessary. An assistant is required for the adjustment work. Due to the large number of setting options for the standing device, some of the settings described may not be possible for some users or under certain circumstances.

Check that all connections are tight after adjustment work!

The cushion of the table top must never be set below the top edge of the pool!

6.1 Body size

Adjustment to the patient's height is made via the extension height of the table tubes. The correct setting is selected when the therapy table is at arm height or slightly higher if required for a standing patient. The table height adjustment is supported by two built-in gas springs.

Loosen the clamping levers on the standpipe and pull it out as required (see Figure 3).

The clamping lever must always be tightened when the standing device is in use!

Figure 2: Adjusting the table height


Figure 3: Clamping lever for table height adjustment

6.2 Back pad or hip pad

An assistant is required to use the back cushions!

The height of the back cushion or hip support is adjusted by loosening the clamping levers on the main tubes and the depth is adjusted by loosening the clamping levers on the rear cross tube. Adjustment to the width is made using the catches on the rear cross tube. The pads are positioned above the buttocks in the area of the lumbar vertebrae when the patient is standing (see Figure 4). They should have a stabilizing effect in conjunction with the heel support and knee pad.

The clamping lever must always be tightened when the standing device is in use!

Figure 4: Position of the hip pad in the buttocks

6.3 Foot and heel holder

Adjust the heel holders by repositioning them on the perforated grid plate so that the user has a firm hold (4 positions are possible).

The patient can get in by opening the Velcro straps on the harness. The straps are closed tightly over the instep. The Velcro straps serve as toe fixation. (see Figure 5)

Check that the foot is firmly attached!


Figure 5: Adjusting the heel switches and toe fixation

6.4 Knee pads

The height and depth of the knee pads are adjusted by loosening the vertical and horizontal tube catches. The knee pads are correctly adjusted when they are positioned below the knee on the shin. The knee pads should first be preset and corrected if necessary when the patient is standing. They stabilize the patient's stance in conjunction with the heel support and back pad. (see Figure 6)

Figure 6: Adjusting the knee pads

6.5 Upper body pad

(Optional as accessory item number KR-630-X)

6.5.1 Adjusting the backrest cushion

The back pad supports the middle of the back at shoulder blade level. To attach the pad, turn the knob on the holder a full turn and pull it out. Slide the pad holder into the bracket with the holes facing outwards. Lower the tube into the holder and release the knob as soon as the breast pad is at the desired height. Let go of the button so that the integrated

Indexing pin snaps into the nearest hole in the tube. By tightening the knob, the backrest cushion is fixed in position.

6.5.2 Lateral armrest

The chest pad is supplemented by side pads. Two width-adjustable pads can be attached to the chest pad, which are positioned on both sides of the chest at the height of the armpits. The height of the pads depends on the height adjustment of the chest pad and should be adjusted so that there is a distance of three fingers between the top of the side chest pads and the armpit. See previous instructions for adjusting the height of the breast pads.

6.5.3 Adjusting the width of the lateral hip supports

The hip supports can be individually adjusted in width to center the thorax. Their position is held with locking pins. Pull out the plungers and slide the arm supports along the horizontal tube to the desired position. Snap the plungers into the nearest holes and make sure that they are locked in the desired position. Ensure that the supports are positioned symmetrically.

6.5.4 Attaching and adjusting the headrest

The headrest supports the head at the level of the back of the head. To attach the headrest, pull out the locking pin on the chest support. Slide the headrest tube into the holder with the holes facing outwards. When you have reached the desired height, release the indexing plunger so that it engages in the nearest hole in the headrest tube.

The depth of the headrest should be adjusted so that the head is supported but not pushed forward. To adjust the depth of the headrest, pull the plunger on the top of the headrest upwards and slide the headrest in or out to the desired position. When you have reached the desired position, release the plunger and make sure that it engages in one of the holes in the headrest.

Figure 7: Upper body pad

6.6 Breast pad

The chest support supports the patient when standing. After centering and aligning the chest support on the table, it can be attached to the table of the standing device using two screw clamps. (see Figure 8)


Figure 8: Breast pad

6.7 General instructions for use when standing up with the Levatrain

Put on the patient harness before approaching the standing device.

The patient belt should be positioned so that its lower edge is clamped between the buttocks and the wheelchair seat. Then fasten the abdominal belt (see Figure 9).

Swivel the legrests or footrest to the side. Place the feet in the foot guides. The heel sections can be removed for this purpose. Drive the wheelchair so close that the front wheels touch the recesses provided in the footplate. Secure the wheelchair against rolling away. (see Figure 5)

Figure 9: Patient restraint belt

Once the control unit has been connected to the mains using the mains cable, the device is ready for operation. The drive is automatically switched off in the end position of the two belt tongue pull cords.

The additional forced shutdown is triggered as long as the black emergency stop button is held down. The appliance can only be operated again once it has been released. The push button is located underneath the tabletop. (see Figure 10)


Figure 10: Position of emergency stop pushbutton

Each upward or downward movement of the straps is effected by pressing the respective button on the handset. The buttons are marked with arrows according to their function. (see Figure 11)

Figure 11: Manual switch for Levatrain and Levatrain Max

To give the wheelchair more room to move, the castors can be adjusted outwards using a screw joint. (see Figure 12)


Figure 12: Adjusting the castors

6.8 General instructions for use when standing up with the Levaflex

The standing device with tilting function (Levaflex) offers you more comfort when getting out of bed than the standard version. The armrests, knee support and footplate are tilted towards the patient before standing up. An additional tilt motor is used for tilting.


Figure 13: Hand control for Levaflex

Figure 14: Dismantling process

Procedure for getting up		
	Description	Button on the handset
1	Tilt the standpipes back by pressing the hand control. Click in both pull ropes with the patient restraint belt in place.	
2	Tension the traction ropes by pressing the button until the ropes are taut but the patient is not yet lifted.	
3	Pressing the stand-up function button activates both motors simultaneously and starts the stand-up process. (see figure 14)	

Table 3: Procedure for standing up

Reset procedure:		
	Description	Button on the handset
1	If you press the button on the hand control, the lowering process is carried out, both drive units move simultaneously and the patient is gently brought into the sitting position.	
2	Press this button to release the tension on the pull cords	
3	If the patient with the wheelchair is removed from the standing device, the table can be brought into the horizontal position by pressing this button.	

Table 4: Reset procedure

6.9 Behavior in the event of a power failure

If the power fails during use of the electric belt retractor, the user must call an assistant to help the patient back into the wheelchair.

The following description assumes that the assistant is right-handed. For a left-handed person, proceed in reverse.

First, the wheelchair is turned around the left front wheel just enough so that the assistant can stand behind the patient, but it is still possible to place the patient in the wheelchair. The assistant stands behind the user, opens the Velcro fastener, grasps the patient's chest with the right arm under the patient's right armpit and opens the left buckle of the patient restraint belt with the left hand (if necessary, relieve the patient restraint belt by applying some pressure via the body buckle). The assistant's left hand now grasps under the patient's buttocks to support them as they slide back into the wheelchair.

6.10 Operating the brake rollers

The standing device has four swivel castors with individual brakes. The standing device should always be braked at the installation site and during use using the castor lock. (see Figure 15)


Figure 15: Exemplary rollers in braked and unbraked condition

Standing training with unbraked castors can lead to injuries to the patient.

7 Ambient conditions

7.1 Storage conditions



Storage temperature	min. -10 °C max. +50 °C	
Relative humidity	min. 20 % max. 80 %	
Air pressure (at altitude ≤ 3000 m)	min. 700 hPa max. 1060 hPa	

Table 5: Storage conditions

7.2 Operating conditions




Operating temperature	min. +5 °C max. +40 °C	
Relative humidity	min. 20 % max. 80 %	
Air pressure (at altitude ≤ 3000 m)	min. 700 hPa max. 1060 hPa	

Table 6: Operating conditions

8 Technical data subject to change!

Designation	KR-204-B/S	KR-214-B/S	KR-216-B
Nominal voltage	~230 V/50Hz		
Rated power	690 VA		
Device type B according to IEC 601-1			
Protection class			
Sound power level	63 dB(A)		
IP protection class Hand switch, control box, Motor	IPX4		
IP protection class Standing device lifting mechanism	IPX1		
Duty cycle ED 10%	maximum 6 minutes/hour		
max. patient weight	120 kg	120 kg	200 kg
Patient size	140 - 200 cm	140 - 200 cm	140 - 200 cm

Empty weight	approx. 62 kg	approx. 65 kg	approx. 65 kg
Width	78 cm	78 cm	89 cm
Length	100 cm	100 cm	110 cm
Table height	100-132 cm	100-132 cm	100-132 cm

Table 7: Technical data

Repairs may only be carried out by ISKO specialist personnel or by persons authorized and trained by ISKO with comprehensive product knowledge. Any warranty and liability claims will be rejected if this provision is not complied with.

9 Materials used

The upright unit is manufactured as a welded tubular steel construction. The surfaces are powder-coated or galvanized. All wooden parts are either laminated or lacquered. The surfaces of this product are harmless to the skin from a health point of view.

10 Service and care

All household cleaners without ammonia and abrasive cleaners are permitted for cleaning the tube parts, the lying surface and the wooden parts with a damp cloth. Solvents (e.g. nitro) destroy the coating of the tubes and the lying surface!

Mechanical cleaning (e.g. scraping, sanding) or blast cleaning of the standing device is not permitted. All pivot points of the moving parts, including the bearing eyes on the adjuster, are fitted with maintenance-free plain bearings and must not be oiled or greased.

11 Useful life of the product

With an expected average level of use in home care, the service life of the standing device is 10 years. Lack of maintenance and excessive strain on the product can significantly reduce the service life of the standing device. The expected service life in professional nursing home use is 7 years.

12 Disinfection

- To ensure that the standing device functions perfectly, every ISKO standing device should be cleaned, disinfected and checked after each use so that it can be used again immediately.
- Improper cleaning/disinfection of the standing appliance can be dangerous.
- Depending on the degree of soiling, we recommend cleaning the standing appliance with a damp cloth or similar.
- For the wiping and spray disinfection are disinfectants in their intended concentration. (see manufacturer's instructions)
- The dilution ratio recommended by the manufacturer in the respective instructions for use must be used.

Solvents are not permitted.

Abrasives or scouring pads must not be used.

12.1 Specifications of the cleaning agents and disinfectants:

- The working solutions should normally be used freshly prepared.
- Do not exceed or fall below the specified concentrations.
- They must not contain any corrosive or caustic components.
- They must not contain any substances that change the surface structure or the adhesion properties of the materials.
- Lubricants must not be attacked by the cleaning agent and disinfectant.

Under no circumstances should soap or detergents be added to the disinfectant. There is a risk of explosion and fire if alcohol-based products are used over large areas.

The use of unsuitable detergents and disinfectants may cause damage to the surface coating, for which ISKO KOCH GmbH cannot be held liable.

13 Malfunctions and their rectification

	Error	Measure
1	None of the motors react to the manual button actuations.	Check the plug connection of the mains connection cable. Check the key switch (ON/OFF)! Check the plug connection between the hand control and motor cables with the control unit.
2	One motor does not respond to the switch actuation	Check plug connection between motor cable and control unit.

Table 8: Malfunctions and their rectification

Modifications, readjustments and repairs to the standing device that cannot be carried out in accordance with the above instructions may only be carried out by the manufacturer directly or by a workshop authorized by the manufacturer.

14 Recommended accessories

Article description	Order number
Support bar with hip supports	KR-620-X
Straight upper body support with head and arm support	KR-631-X
Breast pad 15 cm, is attached to the table with screw clamps	KR-744-X
Breast pad 25 cm, is attached to the table with screw clamps	KR-745-X
Table top horizontally adjustable	KR-710-X
Foot plate horizontally adjustable	KR-720-X

Table edge	KR-602-X
Patient restraint belt, size S, standard	KR-625-X
Patient restraint belt, size M, standard	KR-624-X
Patient restraint belt, size L, standard	KR-623-X
Patient restraint belt, size S, with base	KR-608-X
Patient restraint belt, size M, with base	KR-607-X
Patient restraint belt, size L, with base	KR-606-X
Patient restraint belt, size XL, with base	KR-609-X

Table 9: Recommended accessories

15 Maintenance

15.1 Legal basis

The Medical Device Regulation (EU) 2017/745 (MDR) and national laws and regulations require operators of medical devices to ensure that the medical device is in a safe operating condition throughout its entire service life.

15.2 Maintenance intervals

As a requirement of the Medical Devices Operator Ordinance §4 (Maintenance), a thorough visual inspection (1), a functional test (2) and a functional test (3) must be carried out after the medical device has been in operation for two years at the latest.

(2) and a current leakage test (3) in accordance with DIN EN 62353:2015-10.

(1) Special attention must be paid to the following points during the visual inspection:

- Tight fit of all screw connections
- Mobility of the pivot points
- Check the mains supply cable for crushing or shearing points
- Checking the strain relief of the mains supply cable

(2) Special attention must be paid to the following points during the functional test:

- Function of all electrically operated movements
- Fully extend and retract all motors on the standing device (without patient) until they switch off automatically. (Limit switches in the motors must switch off with an audible click)
- Functionality of the brakes
- Mobility of the triggers
- Checking the manual switch

Functional tests and current leakage tests may only be carried out by ISKO specialist personnel or by persons authorized and trained by ISKO with comprehensive product knowledge.

15.3 Spare parts

All spare parts for this medical device are available from ISKO KOCH GmbH, quoting the serial number, order number and article number (these can be found on the type plate attached to the medical device).

To ensure that functional reliability and any warranty claims remain valid, only original ISKO KOCH GmbH parts should be used for spare parts.

ISKO KOCH GmbH

Egerländer Street 28

95448 Bayreuth

Phone: +49(0)921/150845-0 (Monday - Thursday 8:00 - 17:00 & Friday 8:00 - 16:00) Fax:

+49(0)921/150845-45

E-mail: info@isko-koch.de

15.4 Notes on documentation

According to the Medical Device Operator Ordinance and Medical Device Regulation (EU) 2017/745 (MDR), there is a documentation obligation for:

- Maintenance
- Incidents / near misses

If unusual hazards for the product are foreseeable at the installation site of the standing appliance (supply line is on the ground; children playing; pets; ...), the electrical cables in particular must be constantly checked and suitable measures taken to avert hazards.

16 Re-use

Before using the standing device again, a thorough visual and functional check of all electrically operated functions as well as a current leakage test in accordance with DIN EN 62353:2015-10 must be carried out as described in the **maintenance intervals** section. The points on **service and care & maintenance** specified in the operating instructions must be observed when cleaning the standing device.

17 Waste disposal

17.1 Disposal of the device

The appliance and accessories, if present, should be disposed of in an environmentally friendly manner and in accordance with legal regulations. Please comply with the applicable **waste separation regulations!** If you are still unsure about this, please contact your local authority or waste disposal company.

17.2 Disposal of electrical components

*if electrical components are included in the medical device

In accordance with Directive 2012/19/EU - WEEE2, this medical device is classified as electrical equipment. All electrical components are free of non-authorized substances classified as harmful in accordance with RoHS II Directive 2011/65/EU. In addition, replaced electrical components must be disposed of in accordance with the European directives (see Directive 2012/19/EU - WEEE2).

17.3 Disposal of the packaging

The EU Waste Framework Directive 2008/98/EC is decisive for the handling of packaging disposal. Reusable materials must be recycled in accordance with national regulations.

18 EC Declaration of Conformity

As the manufacturer, we declare under our sole responsibility that our standing appliances comply with the essential requirements of the EC Directive for

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